### PH-MESSSTICK PHT-02 ATC

BEST.-NR. 101137

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ITEM NO. 101137

### INTENDED USE

The probe allows measurement of the pH of aqueous solutions in the range of 0 to 14. The pH value is a measure of the acidity or alkalinity of an aqueous solution. The measurement is applicable to liquids that are not subject to any voltage. Automatic temperature compensation ensures reliable readings even at different measuring temperatures.

The device can be calibrated for the measurement to compensate for wear of the electrode. To this end, the device houses two calibration controls that allow adjustment of the zero point (pH7) and the lower reference point (pH4) or the upper reference point (pH10), if necessary.

The device conforms to protection class IP57 and thus, dust and waterproof, and can be immersed up to max. 1 m for a short time.

Operated via two Lithium button batteries of the CR2032 type. No other electrical power supply is permitted. The measuring stick turns off automatically after 5 minutes.

Measurement is not allowed in areas prone to the risk of explosions (Ex).

Use other than that described above will lead to damage to the product. No part of the product may be modified or converted! The safety instructions are to be observed without fail!

### PACKAGE CONTENTS

PHT-02 ATC pH meter	pH7 buffer solution
Batteries (2x CR2032)	pH4 buffer solution
Alignment screwdriver	plastic box
Soaking solution (3mol KCL)	operating instructions

### SAFETY INSTRUCTIONS AND HAZARD WARNINGS

Please read all the instructions before using this device, they include important information on its correct operation.

Damage due to failure to follow these operating instructions will void the warranty! We do not assume any liability for any resulting damage!

We do not assume any liability for material and personal damage caused by improper use or non-compliance with the safety instructions! In such cases, the warranty will be null and void.

An exclamation mark in a triangle indicates important information in these operating instructions which is to be followed strictly.

→ The "arrow" symbol is used where special tips and notes on operation are provided.

This device has been CE-tested and it therefore meets the necessary national and European guidelines.

### Persons/product

The unauthorized conversion and/or modification of the product is inadmissible for reasons of safety and approval (CE).

The device must be set up correctly. Therefore, please follow the operating instructions carefully.

Measurement instruments used at schools, training facilities, do-it-yourself and hobby workshops should not be handled unless supervised by trained, responsible personnel.

Be careful when handling flammable, caustic or alkaline liquids. In that case, wear safety gloves, goggles and an apron. Only measure in well-ventilated environments.

Do not switch the measuring instrument on immediately after it has been taken from a cold in to a warm environment. The condensation that forms might destroy your device. Allow the device to reach ambient temperature before switching it on.

Do not leave packing materials unattended. It could be a dangerous plaything for children.

Do not expose the device to extreme temperatures, strong vibrations or strong mechanical stress. The correct operation of the display is only guaranteed in the temperature range 0 to +50 °C.

Apart from changing battery or electrode, the housing may not be opened or dismantled.

The glass electrode is well-protected within the replaceable electrode and may not be touched.

Electronic devices are not toys and must be kept out of the reach of children.

Do not leave the buffer solutions unattended. They could be swallowed by children or pets. If swallowed, ask the person to drink plenty of water. The water taken in may not be vomited. Then consult a doctor immediately.

Avoid direct contact of buffer solutions with the skin or the eyes. In case of contact, wash the affected skin area with soap and a lot of cold water. In case of contact with the eyes, wash the eyes under running water for about 15 minutes.

### Batteries

Make sure that the polarity is correct when inserting the batteries.

Remove the batteries, if you are not going to use the device for a while, to prevent damage from leaking. Leaky or damaged batteries can cause burns upon contact with the skin. Therefore, you should wear protective gloves when handling damaged batteries.


Keep them out of the reach of children. Do not leave batteries lying around. There is a risk that the batteries may be swallowed by children or pets. If they are swallowed, consult a doctor immediately.

Replace all of the batteries at the same time. Mixing old and new batteries may lead to battery leakage and thus to damage to the device.

Do not dismantle batteries, do not short-circuit them or throw them into a fire. Never try to charge regular batteries. There is a risk of explosion!

### DESCRIPTION OF THE PARTS

- 1 Screw lid with battery compartment and calibration controls
- 2 Display
- 3 Power button to switch on and off
- 4 Screw fitting for replaceable electrode
- 5 Replaceable electrode
- 6 Glass electrode
- 7 Moistening cap
- 8 Protective cap
- 9 Battery compartment
- 10 Calibration controls for reference point (SPAN)
- 11 Holder for alignment screwdriver
- 12 Calibration controls for zero point (ZERO)

### INSERTING/CHANGING BATTERIES

Insert the included batteries before putting the device into operation for the first time.

A battery change becomes necessary when the display becomes hard to see or when the device can no longer be turned on. When the batteries are empty, the measuring stick turns off automatically within 3 seconds after being switched on.

To insert/replace the batteries proceed as follows:

Screw the screw lid (1) of the device.


Remove the used batteries from the battery compartment (9), if necessary. The opening at the back of the battery compartment helps to lift off the batteries with an object.

Place two identical new batteries with the correct polarity into the battery compartment. The positive pole of the batteries has to point downwards to the device.

Close the battery compartment carefully.

All screw-fitted parts are protected from infiltration of water using rubber seals. Make sure that the seals are not pinched or damaged. Regular use of silicone grease is recommended.

### COMMISSIONING

#### Preparation

The meter is used to measure a glass electrode with an open-pore membrane. This membrane may not dry out, so as to preserve its special properties.

Always make sure that there is sufficient soaking solution in the moistening cap (7). A few drops of soaking solution are sufficient here.

Remove the protective cap (8) from the device. This cap is merely plugged.

Screw the moistening cap (7) of the replaceable electrode (5) and check the liquid. If necessary, fill a few drops of the soaking solution into the moistening cap and screw it again onto the electrode.

Always screw the moistening cap onto the electrode to prevent drying out during longer measurement downtimes.

Provide distilled or deionised water and soft paper towels to clean the electrode.

In order to ensure the function of this device over a long period of time, the pH electrode needs to be cleaned with distilled water after each measurement (see maintenance and cleaning).

#### ALIGNMENT / CALIBRATION

Alignment / calibration should always be performed after an electrode change or before the first measurement. In a measurement series, it is sufficient to perform a calibration after about every 10<sup>th</sup> measurement or after several days of measurement downtime. If the supplied buffer solutions are used up, these can be purchased optionally.

Calibration of the zero and the reference point should always be performed sequentially.

Screw the battery compartment (1) lid of the device.

There is a holder inside with a small alignment screwdriver (12). Remove this. Leave the battery compartment open.

Remove the protective cap (8) and the moistening cap (7).

Rinse the electrode head and the glass electrode with distilled water and pat the replaceable electrode dry.

Do not touch the glass electrode inside the electrode head, as this may damage it.

## Zero point calibration

Take some pH7 buffer solution in a small container. The glass electrode should be completely covered during immersion.

Turn on the meter by pressing the power button (3).

Wait a while until the display has stabilised and then dip the electrode head in the buffer solution. Stir slowly with the meter in the buffer solution until the reading has stabilised. Leave the meter immersed and set the reading to a pH 7.00 using a calibration control "ZERO" (12). To do this, use the small alignment screwdriver.

→ The pH changes slightly at a temperature deviation of 25 °C. The corresponding pH value is specified on the buffer solution container for the different temperatures. Use the corresponding value for alignment.

Rinse the electrode with distilled water, rinse and pat the replaceable electrode dry.

## Reference point calibration

The reference point shall always be close to the measured pH of the liquid. Select the reference point pH4 if you want to measure acidic liquids (< pH7), and the reference point pH10 for alkaline liquids (> pH7). The reference buffer solution pH4 is supplied with the meter.

Take some pH4 or pH10 buffer solution in a small container. The glass electrode should be completely covered during immersion.

Wait a while until the display has stabilised and then dip the electrode head in the buffer solution. Stir slowly with the meter in the buffer solution until the reading has stabilised. Leave the meter immersed and set the reading to a pH 4.00 (or pH 10.00) using the calibration control "SPAN" (10). To do this, use the small alignment screwdriver.

→ The pH changes slightly at a temperature deviation of 25 °C. The corresponding pH value is specified on the buffer solution container for the different temperatures. Use the corresponding value for alignment.

Rinse the electrode with distilled water, rinse and pat the replaceable electrode dry.

Place the alignment screwdriver in the holder (11) within the battery compartment and close the battery compartment carefully.

## MEASURING

Perform a calibration if the meter is not used for a long time or a electrode was changed.

Proceed as follows for measuring:

Pull the protective cap (8) off the device. This cap is merely plugged.

Screw the moistening cap (7) of the replaceable electrode (5) and check the liquid. The electrode must be moist.

Clean the electrode with distilled water and tap it dry. Never touch the glass electrode (6) while doing so as this may damage it.

Turn off the device by pressing the power button (3).

Dip the electrode into the liquid to be measured. The glass electrode has to be completely covered by the liquid.

Stir slowly with the meter in the buffer solution until a stable reading is displayed. Read this from the display.

Turn off the device after the measurement by pressing the power button (3).

Clean the electrode and screw the moistening cap onto the device. If necessary, add a little soaking solution into the moistening cap.

Replace the protective cap (8) on the meter and stow it away carefully in the plastic box.

→ Measurements in hot liquids cause the electrode to wear much faster.

## CHANGING THE ELECTRODE

The electrode is a wear part. Worn electrodes are not covered by the warranty/guarantee.

If measurement is not possible or calibration cannot be carried out, the electrode is usually worn and must be replaced.

To replace, proceed as follows:

Make sure that the meter is switched off and the battery compartment is closed.

Remove the protective cap (8) from the device.

Loosen the screw fitting (4) of the electrode and take it off.

Carefully pull the replaceable electrode (5) off the meter without turning the electrode.

Insert the new replaceable electrode into the meter. Please note the notch in the plug region. This fits correctly on the meter only on the right side.


Screw the replaceable electrode again with the screw fitting (4).

Check the liquid in the moistening cap and refill the soaking solution, if necessary.

Turn on the meter and check the operational readiness. A reading must be displayed. If no reading is displayed, please check the plug connection. Turn the device off.

Perform a calibration before the first measurement.

Replace the protective cap (8) on the meter and stow it away carefully in the plastic box.

## TROUBLESHOOTING

In purchasing the measuring instrument, you have acquired a product with state-of-the-art design that is operationally reliable. Nevertheless, problems or faults may occur. Which is why we would like to describe here how you can solve any problems:

Always follow the safety instructions!

Error	Possible cause	Remedy
The meter cannot be switched on	The batteries are flat.	Change the batteries.
Incorrect readings	The temperature of the liquid is out of the allowable range of 0 to +50 °C.	Bring the liquid within the given temperature range.
	The device has not been calibrated.	Perform a calibration.
Calibration is not possible	The electrode is dirty.	Clean the electrode.
	The buffer solutions are used.	Use new buffer solutions.
	The electrode has dried up.	Recover the electrode.
	The electrode is used up.	Replace the electrode.

Repairs other than described above are carried out exclusively by authorised specialist staff. If you have questions concerning the use of the meter, please contact our technical support team.

## CLEANING AND MAINTENANCE

Maintenance is limited to changing of batteries and electrodes.

Use a clean, lint-free, antistatic and slightly damp cloth to clean the device. Do not use any abrasive or chemical agents or detergents containing solvents.

The electrode must be cleaned after each measurement. Distilled water is sufficient for normal cleaning. A pepsin-cleaning fluid is required for measurements in protein-containing liquids (e.g. milk). Leave the electrode in the pepsin-cleaning fluid for about 5 minutes and then rinse thoroughly with distilled water.

Pat the electrode dry and screw the moistening cap onto the electrode.

All screw-fitted parts are protected from infiltration of water using rubber seals. Make sure that the seals are not pinched or damaged. Regular use of silicone grease is recommended.

## Recovery of the glass electrode

A dried glass electrode has a negative effect on the service life. Dried electrodes may be recovered to render the open-pore glass layer permeable again.

Dip the dried electrode in plenty of soaking solution (3 mol KCl solution) for min. 24 hours and then try re-calibrate. If this is not successful, extend the time again by 24 hours. Then, a measurement should be possible; if not, the electrode must be replaced.

## DISPOSAL

Electronic products are raw material and do not belong in the household waste. When the device has reached the end of its service life, please dispose of it, according to the current statutory requirements, at your local collecting site. It is prohibited to dispose of the device with domestic waste.

## Disposal of used batteries!

As the end user, you are required by law (Battery Ordinance) to return all used batteries/rechargeable batteries; disposal of them in the household waste is prohibited!

Batteries/rechargeable batteries containing harmful substances are marked with the following symbols, which indicate that disposal in the household waste is prohibited. The symbols of the critical heavy metals are: Cd = Cadmium, Hg = Mercury, Pb = Lead. You can return your used batteries/rechargeable batteries free of charge at the official collection points in your community, in our stores, or at places where batteries or rechargeable batteries are sold!

You thereby fulfil your statutory obligations and contribute to the protection of the environment!

## TECHNICAL DATA

Display.....	LCD
Measuring range.....	0.01 - 14.00 pH
Accuracy.....	±0.01 + 1 count
Temperature compensation.....	automatic, 0 to +50 °C
Power supply.....	2 Lithium batteries (CR2032)
Weight.....	approx. 85 g
Dimensions (H x Ø) in mm.....	168 x 37
Weight.....	approx. 85 g

## Legal notice

This is a publication by Conrad Electronic SE, Klaus-Conrad-Str. 1, D-92240 Hirschau (www.conrad.com).

All rights including translation reserved. Reproduction by any method, e.g. photocopy, microfilming, or the capture in electronic data processing systems require the prior written approval by the editor. Reprinting, also in part, is prohibited. This publication represent the technical status at the time of printing.

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© Copyright 2016 by Conrad Electronic SE.