BRESSER®

**Bedienungsanleitung
Instruction Manual
Mode d'emploi
Manual de instrucciones**

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Kaufdatum:.....

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CAUTION!

To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

All parts

- 1 Eyepiece
- 2 Binocular head
- 3 Eye relief control
- 4 Diopter setting
- 5 Microscope stand
- 6 Nosepiece with 4 objectives
- 7 Microscope table (mechanical desk with vernier scale)
- 8 Lengthwise knob
- 9 Cross knob
- 10 Condenser
- 11 Condenser height adjustment
- 12 Coarse focusing wheel
- 13 Fine focusing wheel
- 14 Illumination – On/Off switch
- 15 Dimmer
- 16 Immersion oil, light bulb and fuse
- 17 Fixing screw for binocular head
- 18 up stop for microscope table

1. General/Location:

Before you set up the microscope, you must choose a suitable location. Firstly, you must make sure that your microscope is on a stable and solid surface. For observations with the electric illumination an electrical power outlet is required.

2. Eyepiece

For observation with this microscope, two eyepieces (1) are required. With the eyepiece and with a chosen objective lens, you can choose the magnification of the microscope. The eyepiece is inserted into the head from above.

TIP: Begin each observation with the lowest magnification. So, the centering and focussing of the object to be viewed is easy.

3. Built in electric illumination

The microscope has a built-in electric illuminator. It works by plugging the connection cable into the electricity supply. It can be switched on or off with the switch (14) in the rear of the microscope. Use the dimmer (15) to change the intensity of the light.

4. Adjust the eye relief

The eye relief varies from person to person. In order to achieve a perfect interaction between your eyes and the eyepieces, the eyepiece distance should be adjusted before the first observation. Look through the eyepieces and adjust the eye relief control (3) until you reach a single, round image. Keep the value on the scale in mind for a faster adjustment the next time.

5. Observation

- You begin each observation with the lowest magnification (eyepiece 10x and objective 4x). So that the object to be viewed is centred and properly focussed.
- The higher the magnification the more light is required for good picture quality.
- For applications with the 100x objective, it is recommended that examinations are done only with oil immersion.

Begin with a simple observation.

Turn the nosepiece (6) onto the smallest magnification and insert the 10x eyepieces. Now, place a prepared slide under the objective on the microscope table (7). Move the preparation on the table with the two adjusting screws (8 and 9). Now, look through the eyepieces (1) and

turn the coarse adjustment (12) carefully until the image becomes visible. The definition is applied by the fine adjustment (13). Make sure that you never overtighten the fine adjustment. If you don't reach a sharp image after the fine adjustment, you may adjust the diopter setting (4) until it matches your eye. For this, turn the diopter setting (4) to O, look through the eyepieces and close your left eye. Focus the object for your right eye using the coarse and fine focussing wheels (12 and 13). Now, close your right eye and look only with the left eye through the left eyepiece. If the image isn't sharp, turn the diopter setting (4) to get a sharp image.


6. The microscope table

With the microscope table (7), you can look at your preparation in a precise position and to the exact millimeter. The object is placed between the clamps on the microscope table. Move the object, with help of the axis-adjustments (8 and 9), directly under the objective. With the built-in vernier at both axes you can now specifically set and shift the object. It can now viewed with different magnifications.

7. Changing the bulb (Fig. 2)

Turn the microscope off using the On/Off switch and pull the plug out of the socket. Wait, until the lamp has cooled down surely. Then lay the instrument carefully on its side. Undo the screw (1) to open the housing. The bulb and its mounting are now visible. Carefully remove the old bulb from its socket (2). Insert a new one carefully. Clean the new one with a little alcohol to remove any fingerprints. Then close the housing again and retighten the screw (1). The microscope is now ready for use again.

8. Changing the fuse (Fig. 2)

A fuse is built in for safety reasons. Should it burn through, replace it as below: Pull the plug and lay the instrument carefully on its side. Carefully, undo the "FUSE" screw (3). Use a small screwdriver for the purpose and press the screw a little in to undo it. Remove the old fuse and insert a new one of the same type (1 A, T). Retighten screw and fuse again. The microscope is now ready for use again.

9. Oil immersion

Directions for observation with oil immersion. Put a drop of immersion oil (16) onto the upper lens of condenser and raise the condenser by using the height adjustment (11) until the oil touches the slide from below. Centre your object, with middle magnification, and with the light illumination on. Move the nosepiece with the coarse focusing wheel up high and screw in the 100x objective. Take a drop of immersion oil and place it on the preparation. Move the objective down carefully with the coarse adjustment, until it touches the immersion oil. Look through the eyepiece and focus with the fine adjustment. Do not forget to clean the objective, the preparation, and the condenset after the observation with some alcohol.

10. Care and maintenance:

- 1. Avoid disassembling the microscope!
Because of being a precise optical instrument, do not disassemble the microscope casually. That may cause serious damage to its performance. Do not disassemble the objectives.
- 2. Avoid contamination!
You should make sure that dust or moisture do not come in contact with your microscope. Avoid putting fingerprints on any optical surfaces! If dust or dirt nevertheless have gotten onto your microscope, or the accessories, remove it as described below.
- 3. Cleaning non-optical parts
Disconnect the plug from mains socket before cleaning non-optical parts of the microscope (e. g. the frame or the stage). Firstly, remove dust or dirt from your microscope, or accessories, using a soft brush. Then clean the soiled area with a soft, lint-free cloth; the cloth may be moistened with a mild detergent solution. Be sure that the microscope is dry before using. For cleaning the painted parts or plastic components, do NOT use organic solvents such as alcohol, ether, acetone, xylene or other thinner etc.!
- 4. Cleaning optical parts
In order to enhance the optical quality, the lenses of eyepieces or objectives are coated. They should not be wiped because dry dirt or dust may scratch the coating. It is best to remove the parts to be cleaned from the frame prior to cleaning. Always blow loose dust away first. Use lens tissue of good quality or soft fabric moistened with a small amount of pure alcohol (available in the chemist's shop); wipe the lens surface clean.
- 5. Cleaning the oil immersion lens
After working with oil immersion, the immersion oil should be removed from the lens of the 100x objective at the end of each workday. In order to do this, use a lens tissue moistened with a small amount of pure alcohol (available in the chemist's shop). Don't forget to clean the preparation, too.
- 6. Being not used
After usage, cover the microscope with the dust cap (made of plastic) and place it where it is dry and mouldless. In the case of non-use for a longer time, you should stow the microscope and the accessories in their correct containers again. We suggest the storage of all objectives and eyepieces in a closed container with drying agent.

Remember: A well maintained microscope will keep its optical quality for years and thus maintain its value.

11. Technical data:

Binocular head
 Mechanical desk (with vernier alignment)
 Eyepieces: DIN WF 10x
 Objectives: DIN 4x / 10x / 40x / 100x
 Magnification: 40x -1000x
 Illumination: 6 V / 20 W
 Accessories: Immersion oil, light bulb, fuse

12. Aperture:

The numerical aperture (the unit of measurement for the resolution factor of the objectives) is engraved beside the magnification sign on the objectives.

13. The magnification for microscopes

It is calculated by the formula:
 Eyepiece magnification x objective magnification = total magnification
 Example: 10x-eyepiece x 100x-lens = 1000x total magnification

14. EEC conformity explanation

Meade Instruments Europe GmbH & Co. KG, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

73/23/EEC · 93/68/EEC
Registration No. 15001933001

Product description: Transmission-type microscope
Model: BRESSER Researcher Bino

Rhede, 4th June 2005
Meade Instruments Europe

Helmut Ebbert
Managing director

15. Warranty

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Microscopes which become defective during the warranty period can be returned to the dealer where the microscope was bought. The repaired microscopes or a new pair will then be returned to you, free of charge. In the case of defects which occur after the end of the warranty period, the microscopes can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:
Make sure to return the microscopes carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This guarantee does not imply any restriction of your statutory rights.

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