
BRESSER


MicroSet 50x-1200x

Art. No. 88-51000

- DE** Bedienungsanleitung
- GB** Operating Instructions
- FR** Mode d'emploi
- NL** Handleiding
- IT** Istruzioni per l'uso
- ES** Instrucciones de uso
- PT** Manual de utilização

Dear parents,

This product is ideal for children wanting to explore their world in a completely new way. The device is as such, easy to use and care for, rugged and good-looking.

More important to you and of course to us is that it is safe to use. During manufacture, we made sure that this product is as safe it can be for children to use. Some residual risk is, however, unavoidable. This product, after all, is not a toy in the usual sense but rather an optical instrument that children can use to experiment, research and discover their world.

That's why we request your cooperation here. These operating instructions were written for children but please read them through together with your child or children and answer his/her/their questions. Don't forget to explain possible risks.

These are summarised under the heading „warnings“. Please adjust/set up the device together with your child or children and never allow any child to use any of our optical products unsupervised.

We hope all users and their parents will enjoy our products.

Your Bresser team

Dear junior researcher,

Congratulations on becoming the proud owner of this product.

You'll be amazed when reading these instructions just how much can be done and explored with your new device.

Take a look and emerge yourself into the adventurous world of nature and discovery.

It really is exciting and a lot of fun discovering the world with this product.

Before you get started, read the operating instructions fully, as there are a few things you need to know to get the best out of your new device.

The „Warnings“ should be read carefully. Use the product exactly as per the operating instructions to avoid any risk or injury. Keep these instructions in a safe place for later reference. If you give the device away or make a present of it make sure these instructions accompany it.

And now it just remains to say, „Have loads of fun researching and discovering“

Pia

RISK to your child.

Aids with sharp edges and tips are sometimes used with this device. Please store the device and all of its accessories and aids out of the reach of children. There is a risk of INJURY.

This device contains electronic components that are powered by either a mains connection or batteries. Never leave a child unsupervised with this device. The device should only be used as per these instructions otherwise there is a serious RISK of ELECTRICAL SHOCK.

Batteries should be kept out of children's reach. When inserting batteries please ensure the polarity is correct. Leaking or damaged batteries can cause injury if they come into contact with the skin. If you need to handle such batteries please wear suitable safety gloves.

Children should only use this device under supervision. Keep packaging materials (plastic bags, rubber bands, etc.) away from children. There is a risk of SUFFOCATION.

The chemicals and liquids provided should be kept out of reach of children. Do not drink the chemicals! Hands should be washed thoroughly under running water after use.

In case of accidental contact with the eyes or mouth rinse with water. Seek medical treatment for ailments arising from contact with the chemical substances and take the chemicals with you to the doctor.

FIRE-/ DANGER OF EXPLOSION!

Do not expose the device to high temperatures. Use only the mains adapter supplied or those battery types recommended. Never short circuit the device or batteries or throw into a fire. Exposure to high temperatures or misuse of the device can lead to short circuits, fire or even explosion!

RISK of material damage.

Never take the device apart. Please consult your dealer if there are any defects. The dealer will contact our service centre and send the device in for repair if needed.

Do not subject the device to temperatures exceeding 60 C.

TIPS on cleaning

Remove the device from its energy source before cleaning (remove plug from socket / remove batteries)

Clean the exterior of device with a dry cloth. Do not use cleaning fluids so as to avoid causing damage to electronic components.

Clean the lens (objective and eyepiece) only with the cloth supplied or some other soft lint-free cloth (e.g. micro-fibre). Do not use excessive pressure - this may scratch the lens.

Protect the device from dust and moisture. Store the device in the bag supplied or in its original packaging. Batteries should be removed from the device if it is not going to be used for a long period of time.

DISPOSAL

Dispose of the packaging material/s as legally required. Consult the local authority on the matter if necessary.

Do not dispose of electrical equipment in your ordinary refuse. The European guideline 2002/96/EU on Electronic and Electrical Equipment Waste and relevant laws applying to it require such used equipment to be separately collected and recycled in an environment-friendly manner. Empty batteries and accumulators must be disposed of separately. Information on disposing of all such equipment made after 01 June 2006 can be obtained from your local authority.

EEC conformity explanation

 Meade Instruments Europe GmbH & Co KG, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 50088:1996
EN 71 Part 1:1988
EN 71 Part 2:1993

Product description: MicroSet
Model: GMPZ-C900

Rhede, 01-01-2009

Meade Instruments Europe GmbH & Co. KG

Helmut Ebbert
Managing director

All Parts

- 1 Zoom eyepiece
- 2 Focus knob
- 3 Objective turret
- 4 Stage
- 5 Mirror
- 6 Electronic light source
- 7 Base with battery compartment
- 8 Rubber base
- 9 Microscope arm
- 10 Microscope "Biotar"
- 11 Prepared slides
- 12 Cover slips
- 13 Slide
- 14 Petri dish
- 15 Micro projector
- 16 Replacement bulb
- 17 3x magnifying glass
- 18 Settling cylinder
- 19 Reservoir
- 20 Shrimp farming accessories
- 21 Microscope instruments
- 22 Micro-slicer

Your Microscope

What is a Microscope?

The microscope is made of two lens systems: the eyepiece and the objective. We're presenting these systems as one lens each, so that it's easier to understand. In reality,

however, the eyepiece (1) and the objective in the turret (3) are made up of multiple lenses.

The lower lens (objective) magnifies the prepared specimen (e.g. 11), which leads to a magnified image of the specimen. The picture, which you can't see, is magnified once more by the second lens (eyepiece, 1). Then, you can see the "microscope picture."

Assembly and Location

Before you start, choose an ideal location for using your microscope. On the one hand, it's important there is enough light (normal observation). On the other hand, it's important that you can quickly make the room dark when you want to (projector observation). Furthermore, I recommend that you place the microscope on a stable surface, since you won't be able to get any satisfactory results on a shaky surface.

Normal Observation

For normal observation, place the microscope in a bright location (near a window or desk lamp, for example). Remove the microscope from the package and tilt the microscope arm (9) so that it's in a comfortable viewing position for you.

The focus knob (2) is turned to the upper stop and the objective turret (3) is set to the lowest magnification (display: 50x-100x).

Now, take a look through the eyepiece

and adjust the mirror (5) so that you get a uniformly bright circle of light. Or, you can use the lamp. You'll find further tips about the lamp in the next chapter. Now, place a prepared slide (e.g. 11) under the clips on the stage (4), directly under the objective. When you take a look through the eyepiece, you can see the magnified specimen. At this point, it might still be a slightly fuzzy picture. The sharpness is adjusted by slowly turning the focus knob (2). You can now select a higher magnification by turning the objective turret and selecting a different objective.

In doing so, note that the sharpness of the picture must be adjusted again for the higher magnification. Also, the higher the magnification, the more light you will need for good illumination of the picture.

TIP: Avoid positioning the mirror so that it is exposed to direct sunlight, since it can cause a glare, which won't allow you to get a clear picture.

Electric Illumination

To use the electric light source (6), you need two 1.5 V batteries, which are installed in the base of the microscope (7). Remove the rubber base (8) on the microscope and insert the batteries according to the displayed +/-

information. Reattach the rubber base to the microscope.

You can turn on the light by rotating it in the direction of the stage (4). Take a look through the eyepiece and adjust the light so that you get optimal picture brightness. Now you can observe in the same way as described in point 1.3.

The bulbs for the electric light can be replaced. One replacement bulb (16) is included. When you use another bulb, make sure to pay attention to the maximum wattage, which is displayed on the lamp body.

TIP: The higher the magnification used, the more light is required for a good illumination of the picture. Therefore, always start your experiments with one of the lower magnifications.

Projector

The included micro-projector (15) is ideal for looking at a specimen with multiple people at the same time. To use the projector, you absolutely must use the electric illuminator. The best image quality is achieved at the magnifications between 80x and 200x. This way, you can show your results to your friends and family.

Unscrew the zoom eyepiece (1) from the microscope and place the micro-projector on

the open eyepiece support so that the matt screen is pointing behind you. Adjust the light so that the matt screen is lit up brightly.

You can adjust the focus of the prepared slide image using the focus knob (2). Make sure to swivel the microscope arm (9) back into an upright position so that you can comfortably observe the specimen. In a dark room, it is also possible to project the microscope picture on a white wall. To do this, remove the matt screen. Swivel the tripod so that the upper edge of the projector is horizontal. Keep in mind that the image brightness is reduced and the picture will be darker if the projector is far away from the wall.

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