

- **Smartphone Microscope • 2x - 16x**
- **Microscope pour smartphone • 2x - 16x**
- **Microscopio Smartphone • 2x - 16x**
- **Smartphone-Mikroskop • 2x - 16x**
- **Microscópio Smartphone • 2x - 16x**
- **Microscopio per smartphone • 2x - 16x**

Thank you for purchasing a **Smartphone Microscope**. Please read these user instructions carefully before use as they contain important information on the correct use of the **Smartphone Microscope**. Also, keep these instructions for future reference.

Product contents

- Main unit
(with light box and lens plate B already installed)
- Replacement lenses
(Lens plate A and lens plate C)
- Blank slides × 2
- Prepared slide × 1 (includes 3 natural fibres)
- User instructions

Warning:

- Only use this product for its intended purpose. For safety, please read these user instructions carefully.
- Make sure to store this product outside the reach of small children.
- Small children should only use this product under the supervision of a guardian.
- Do not leave the main unit or lens plates in direct sunlight. Lenses can focus the sun's rays and melt the main unit or create a fire hazard.
- Do not attempt to view the sun through the lenses. There is a risk this will cause damage to your eyes.
- This product has small parts that can cause choking. Take care to avoid the accidental swallowing of such parts.
- Take care to avoid water getting onto the main unit or battery box, which can cause the product to malfunction.
- This product is not waterproof. Take appropriate care to avoid spillage when examining liquids.
- Do not disassemble this product. Disassembly can cause the product to malfunction.
- This product uses 2 AAA alkaline batteries. Improper battery use can result in battery overheating, rupture, and leakage.

Please note the following:

- Rechargeable batteries cannot be used as they provide insufficient voltage.
- Ensure batteries are inserted in the correct +/- (plus/minus) orientation.
- Should fluid leaked from a battery enter your eye, immediately rinse the eye with plenty of water and consult a doctor.
- If fluid leaked from a battery comes into contact with your skin or clothing, immediately wash the area.
- Remove the batteries when the product is not in use for an extended period of time. Also, do not use a mixture of new and used batteries.

Main unit specifications


Dimensions: W70 × H130 × D27 (not including protruding parts)

Batteries: 2 AAA alkaline batteries (Do not use rechargeable batteries)

Magnification (when combined with smartphone camera digital zoom)

- Lens plate A: approx. 2 to 10 times
- Lens plate B: approx. 4 to 20 times
- Lens plate C: approx. 16 to 80 times

* Smartphone model: Magnification and zoom range may vary depending on your camera software and screen performance.

Instructions for basic use

For safety, please follow the directions as presented in these user instructions.

1 Insert batteries into the light box

The LED light requires 2 AAA alkaline batteries. Use a screwdriver to loosen the screw in the battery case at the back of the main unit, open the battery cover, and insert the batteries. Confirm the plus (+) and minus (-) battery terminals are in the correct orientation. Rechargeable batteries cannot be used as they provide insufficient voltage.



2 Operate the switch on the light box

Operate the LED light switch on the light box to confirm the light is working. Leave the light box attached to the main unit when examining translucent specimens.

3 Install a lens plate in the main unit

There are 3 different lens plates. The product comes with the medium-magnification lens plate B already installed. Change the installed orientation of the lens plate to match the position of the camera on your smartphone.

4 Use the focus adjustment thumbwheel to raise the main unit slightly

Turn the focus adjustment thumbwheel on the side of the main unit and confirm it raises and lowers the smartphone platform. Then, turn the thumbwheel to raise the smartphone platform slightly.

5 Place the smartphone on the Smartphone Microscope

Activate the camera mode on the smartphone, and hold the smartphone so its camera lens is aligned above the lens on the Smartphone Microscope. Slowly lower your smartphone down onto the platform, making sure to keep the round lens of the Smartphone Microscope in the center of the smartphone display.

6 Place the specimen under the lens

Place the specimen you wish to examine onto a blank slide, and adjust the position of the slide so your specimen is aligned under the Smartphone Microscope lens.

7 Adjust the focus

Turn the focus adjustment thumbwheel to bring the specimen into focus. When the specimen is almost in focus, adjust the brightness and focus settings on the smartphone. If this does not bring your specimen into proper focus, try turning the focus adjustment thumbwheel again.

8 Capture the image with the smartphone

You can use the smartphone display operation to zoom in and further increase the size of the image. Capture either a still image or video with the smartphone. For further details, please refer to the user manual for your smartphone.


Instructions for applied use

9 Change magnification by changing the lens


At purchase, the medium-magnification lens **plate B** is already installed in the main unit. You can enjoy a wider range of magnifications by changing the lens plate. Image quality at the periphery of the image decreases and peripheral distortion increases at higher magnifications.


Depending on the smartphone used, peripheral areas may appear dark and difficult to see. This phenomenon arises due to the wide-angle camera lens used on the smartphone. The images show table salt captured at approximately 10 x, 20 x, and 80 x magnification.

Approx. 10 x


Approx. 20 x

Approx. 80 x

10 The specimen-lens distance that provides the best focus differs depending on the lens

The specimen-lens distance for best focus is about 15 mm from the lens when using lens plate A (image 1), about 8 mm using lens plate B (image 2), and about 2 mm using lens plate C (image 3). When using lens plate A, detach the light box from the main unit. Also detach the light box when examining large specimens that are otherwise too close to the lens for proper focus.

11 Detaching the light box

The light box is hooked onto the main unit. Detach the light box from the main unit by turning the focus adjustment thumbwheel to separate the light box from the main unit then sliding the light box up.

12 Detach the light box to illuminate specimens from the side

When examining an opaque specimen, detach the light box to use it to illuminate the specimen from the side. Alter the position and angle of the illumination, vertical or horizontal, left and right, to best illuminate the specimen for easy examination.

13 Examining large specimens

Detach the light box to enable the examination of large specimens that do not fit under the lens.

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