NUK®

Thermometer Flash

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Congratulations on your purchase of this high-quality NUK product!

The NUK Flash Thermometer will provide you with a quick and reliable way of measuring body temperature on the forehead, as well as a simple way of determining room or bath water and surface temperature for example.

Please carefully read the following instructions for use and keep them in case you want to refer to them later or pass them on to a third party with the device.

Contents:

1	Overview	28
1.1	Use in accordance with purpose.....	28
1.2	Warning and safety instructions.....	29
1.3	Scope of supply and description.....	33
2	Use	33
2.1	Before initial use.....	33
2.2	Basic information.....	34
2.2.1	Fever	34
2.2.2	Fever thermometer.....	34
2.3	Using the thermometer	35
2.3.1	Measurement on the forehead.....	36
2.3.2	Further functions	38
3	Instructions for use.....	41
3.1	Changing the battery.....	41
3.2	Cleaning and storage.....	42
3.3	Disposal.....	42
4	Troubleshooting	43
5	Technical data.....	44
6	Warranty.....	51
7	EU declaration of conformity.....	52

1 OVERVIEW

1.1 USE IN ACCORDANCE WITH PURPOSE

The NUK Flash Thermometer is designed for non-contact measurement of human body temperature on the forehead. Its infrared technology also enables the temperature of surfaces (such as bath water and baby food, for example) to be determined simply and reliably. You need to be aware that this mode is a surface measurement. Select suitable other methods for determining core temperature. Room temperature can also be measured. With the battery supplied, the NUK Flash Thermometer is immediately ready for use. It also has a storage function to enable temperature to be observed over time.

The NUK Flash Thermometer is intended solely for private use indoors, not for commercial use in surgeries or hospitals. The NUK Flash Thermometer was developed for use by adults and young people. Such people must be in possession of full physical, sensory and mental capacity and have sufficient knowledge of how to use infrared thermometers. Children may not use the device, young people only if they have received instruction beforehand.

Under no circumstances is the NUK Flash Thermometer a substitute for examination by a doctor.

Use of the product other than that described here or any modification to the product is prohibited and may lead to injuries to the patient or operator and/or to damage to the product. The manufacturer accepts no liability for damage resulting from use not in accordance with purpose.

Despite use in accordance with purpose and careful handling in accordance with the information in the instructions for use, it is possible that more than 3 years after initial use, your NUK Flash Thermometer may start measuring inaccurately. Faulty temperature results could lead to incorrect diagnosis and therefore have health implications. In this circumstance, do not continue using this device.

Please read the instructions for use in full before using the device for the first time.

Caution

IP22 Degree of protection with regard to water and dust

Symbol for permitted temperature range when storing the device

Representative officer in the European Union

Manufacturer

Type BF medical device

0120

CE mark and code for the appointed office which confirmed conformity with Annex II of Directive 93/42/EEC pertaining to medical devices.

Battery disposal

The device is subject to European Union regulations relating to the disposal of electrical and electronic waste. To protect the environment, the used device should be disposed of at the collection centres provided for it

1.2 WARNING AND SAFETY INSTRUCTIONS

The following provides an overview of all the safety instructions which must be followed to ensure safe use of your NUK Flash Thermometer. In the instructions for use, you will find these safety instructions in the sequence they occur, marked with the symbol  and the number. Notes on effective and optimized use of the device are marked with Note. In these instructions for use, steps to be completed in consecutive order are indicated by a white number in a black circle (e.g. Step 1 as: .

Warning and safety instructions	Section	 no.
Despite use in accordance with purpose and careful handling in accordance with the information in the instructions for use, it is possible that more than 3 years after initial use, your NUK Flash Thermometer may start measuring inaccurately. Faulty measurements of fevers may lead to incorrect conclusions and thus to damage to health. Do not continue to use the device.	1.1	1
Only operate and store the device in accordance with the specified ambient conditions.	2.2.2 3.2	2 20
Only use the NUK Flash Thermometer indoors. External influences such as wind and outdoor temperatures may result in an inaccurate measurement.	2.2.2	3
Check the device for visible damage before using it each time. If you discover any damage, do not use this device.	2.3	4
Check the function and accuracy of the device after any heavy impact (e.g. falling from a height of more than 1 metre). To do this, perform several measurements consecutively. If you have any doubts concerning the results, contact your retailer.	2.3	5
Changes may only be made to the device with the consent of the manufacturer.	2.3	6

Warning and safety instructions	Section	 no.
Use only accessory parts expressly approved for the device.	2.3	7
Keep the protective cap out of the reach of children. There is a risk of choking.	2.3	8
Do not hold the NUK Flash Thermometer in your hands for any longer than necessary, otherwise the device will warm up and may read a body temperature which is too low when checking for fever.	2.3	9
Before each measurement, ensure that the reflective mirror and the sensor lens are clear and clean, otherwise the measurement may be inaccurate.	2.3	10
Use of the device is not a substitute for examination by a doctor. Find out the values which are critical to you and consult a doctor as soon as these values are exceeded. The same applies to measuring results which are ambiguous.	2.3.1	11
The person having their body temperature taken should avoid any fluctuating ambient temperatures for at least 5 minutes before the measurement. You should also avoid any strenuous physical activity or taking a bath for at least 30 minutes before the measurement.	2.3.1	12

Warning and safety instructions	Section	 no.
Make sure that the forehead is dry and free from dirt and cosmetics during the measurement. Avoid placing the device on any scarred areas.	2.3.1	13
If a fever alarm is indicated during the temperature measurement, keep body temperature under continuous observation and, if the person's condition is uncertain, consult a doctor.	2.3.2	14
Measurement in "surface temperature" mode does not give any indication of body temperature and may not be used for diagnosing fever.	2.3.2	15
Do not use the result for surface temperature measurement to assess temperature inside a substance. It may be much higher. Incorrect conclusions may lead to burns or scalds.	2.3.2	16
Keep the battery out of the reach of children. There is a risk of choking.	3.1	17
The battery must be replaced as soon as the "battery empty" symbol appears in the LCD display.	3.1	18
Do not immerse the device in water or expose it to direct moisture. Avoid the reflective mirror steaming up.	3.2	19

Warning and safety instructions	Section	 no.
Do not dispose of used batteries in normal waste. Take them to an in-store battery collection system or to your nearest recycling centre for special refuse.	3.3	21
Do not dispose of the device as normal waste. Instead, take it to an official collection centre for electrical devices.	3.3	22
No modification of this equipment is allowed.	6.	23

1.3 SCOPE OF SUPPLY AND DESCRIPTION

NUK Flash Thermometer as shown in Figures 1 and 2:

- | | | | |
|---|-------------------------|-----------------------------------|-------------------------------------|
| 1 | Protective cap | 6 | Reflective mirror with sensor lens |
| 2 | LCD display | 7 | Unlatching opening for battery |
| 3 | SCAN button | 8 | ON/MEM button |
| 4 | Probe head | 9 | CR2032 battery (inserted in device) |
| 5 | Battery compartment lid | Instructions for use (not shown). | |

2 USE

2.1 BEFORE INITIAL USE

Remove all packaging material from the device, as well as the battery insulation strips. Clean the device as shown in section 3.2.

2.2 BASIC INFORMATION

2.2.1 FEVER

Fever is a state of elevated body temperature.

Body temperature:	"normal" is considered to be around 36.5 to 37.0 °C (97.7 to 98.6 °F). The highest values are measured in the early evening. Measurements differ according to where the temperature is taken (mouth, armpit, anus, ear, temple, forehead).
-------------------	--

Causes:	fever may be due to colds, other illnesses such as respiratory infections, some intestinal infections, hyperthyroidism or the presence of tumours.
---------	--

Symptoms and diagnoses:	sweating, shivering or chills are indications of an elevated body temperature.
-------------------------	--

You should contact a doctor:	when a fever starts to have a debilitating effect on the organism. This depends heavily on the age, constitution and current state of health of the person affected.
------------------------------	--

2.2.2 FEVER THERMOMETER

The NUK Flash Thermometer measures body temperature on the forehead by means of an infrared sensor; no skin contact is necessary.

In forehead mode, the practical fever warning display indicates that an elevated temperature has been measured. In other words, if this reaches ≥ 37.5 °C (99.5 °F), one long and three short consecutive beeps are heard.

Only operate and store the device in accordance with the specified ambient conditions.

Only use the NUK Flash Thermometer indoors. External influences such as wind and outdoor temperatures may render the measurement inaccurate.

The following factors may also have an unfavourable impact on measurement at the forehead:

- skin thickness/skin quality (e.g. scars)
- sweat
- vasoconstrictive medication
- skin irritation
- moisture/iser on the skin

2.3 USING THE THERMOMETER

Check the device for visible damage before each use. If you discover any damage, do not use this device.

Check the function and accuracy of the device after any heavy impact (e.g. falling from a height of more than 1 metre). To do this, perform several measurements consecutively. If you have any doubts concerning the results, contact your retailer.

Changes may only be made to the device with the consent of the manufacturer.

Use only accessory parts expressly approved for the device.

Keep the protective cap out of the reach of children. There is a risk of choking.

Do not hold the NUK Flash Thermometer in your hands for any longer than necessary, otherwise the device will warm up and may read a body temperature which is too low when checking for fever.

Before each measurement, ensure that the reflective mirror and the sensor lens are clear and clean, otherwise the measurement may be inaccurate.

Note

The device must stay in stable ambient (room) temperature for 15 minutes before operating.

2.3.1 MEASUREMENT ON THE FOREHEAD

Use of the device is not a substitute for examination by a doctor. Find out the values which are critical to you and consult a doctor as soon as these values are exceeded. The same applies to measuring results which are ambiguous.

The person having their body temperature taken should avoid any fluctuating ambient temperatures for at least 5 minutes before the measurement. You should also avoid any strenuous physical activity or taking a bath for at least 30 minutes before the measurement.


Make sure that the forehead is dry and free from dirt and cosmetics during the measurement. Avoid placing the device on any scarred areas.

1

Press the ON/MEM button  to switch on the device.

The forehead symbol  will appear on the LCD display. Two beeps will also be heard as signals.

The mode now shown in the display is the standard mode after the device is switched on. The device is now operational.

2

Remove the protective cap from the probe head

Figure 1,
Item 4

3

Keep the device in the centre of the forehead just above the eyes, about 2 to 3 cm away.

Figure 3

4

Briefly press the SCAN button.

A long beep is heard and the egg timer symbol  starts flashing in the display.

As soon as the measurement is complete, this symbol disappears, followed by two short beeps.

5

Read the temperature measured off the LCD display.

If you would like to perform a further measurement, this is possible as soon as the egg timer symbol  has disappeared. To do so, please repeat steps

3**4****5**

, and .

After about 1 minute, if there are no further measurements, the thermometer automatically switches from measuring mode to room temperature mode; this can be seen from the temperature symbol .

Fever alarm

If an elevated temperature of $\geq 37.5\text{ }^{\circ}\text{C}$ ($99.5\text{ }^{\circ}\text{F}$) is measured in forehead mode, one long and three short consecutive beeps are heard.

14


If a fever alarm is indicated during the temperature measurement, keep body temperature under continuous observation and if the person's condition is uncertain, consult a doctor.

Room temperature

If the buttons are not pressed, the device automatically switches to room temperature mode after about 1 minute and now displays room temperature on the LCD display. The temperature symbol  appears.

Mode switch

With the device switched on, press and hold the ON/MEM button  whilst simultaneously briefly pressing the SCAN button (several times in succession). This switches the measuring task between forehead and surface temperature measurement. Either the  or the  symbol will appear on the LCD display.

**Surface
temperature**

In  mode you can determine the surface temperature of objects such as water, milk, clothing or other objects. However, only the radiated temperature of the relevant surface is displayed.

Measurement in "surface temperature" mode does not give any indication of body temperature and may not therefore be used for diagnosing fever.

Do not use the measurement to assess temperature inside a substance. It may be much higher. Incorrect conclusions may lead to burns or scalds.

For measuring in  mode, hold the thermometer as close as possible to the relevant surface.

Note

The further away you move, the larger the range measured and thus also the larger the potential measuring inaccuracies.

If you press the SCAN button briefly, just one measurement is performed and displayed immediately. Press and hold the SCAN button to constantly update the measured temperature and its display.

**Previous
measure-
ment**

This function is valid only with forehead measurement. For every new temperature measurement, the result from the previous measurement will be shown at the top of the LCD display for comparison purposes.

Memory
(25
datasets)


With the device switched on, repeatedly press the ON/MEM button to look at up to 25 previous measurements. The memory is represented by the **(M)** symbol in the LCD display. This function is valid only with forehead measurement.

Switch
between
°C or °F


Changing the temperature display from °C to °F: if the device is in room temperature mode, press and hold the SCAN button whilst simultaneously pressing the ON/MEM button  for 3 seconds.

On the LCD display, the characters switch from "°C" to "°F" and vice versa.

Silent
mode

With the device switched on, press the ON/MEM button  for 3 seconds. The silent mode symbol  starts to flash on the LCD display. Silent mode is activated by releasing the ON/MEM button. Proceed in the same way to deactivate silent mode.

Note With silent mode activated, none of the signal beeps described will be heard.

OFF mode

To switch off the thermometer, press and hold the ON/MEM key  until "OFF" appears in the LCD display. This process may take 8 - 10 seconds.

3 INSTRUCTIONS FOR USE


3.1 CHANGING THE BATTERY

Keep the battery out of the reach of children.
There is a risk of choking.

The battery must be replaced as soon as
this symbol appears in the LCD display:

In this case, proceed as outlined below.

- | | | |
|-------------|--|----------|
| 1 | Take a slim object (e.g. a paper clip) and push it into the unlatching opening of the battery compartment lid whilst simultaneously pushing it outwards with your thumb. | Figure 4 |
| 2 | Firmly holding the device, carefully lift out the battery inside using a slim object (e.g. a paper clip). | Figure 5 |
| 3 | Insert the new CR2032 battery with the positive/+ terminal facing upwards and the negative/- terminal facing downwards by sliding the battery <u>underneath</u> the metal hook and pressing down until it engages. | Figure 6 |
| Note | In the case of incorrect assembly, the battery holder and thus also the device may be damaged. | |
| 4 | Push the battery compartment lid back onto the device until you feel it engage. | |

3.2 CLEANING AND STORAGE

The probe head is the most sensitive part of the NUK Flash Thermometer, so always treat it with particular care. Never touch the reflective mirror or sensor lens inside. To prevent contamination, always store the thermometer with its protective cap on. In order to be able to guarantee accurate temperature measurement, the sensor lens and the reflective mirror must always be clean. If deposits can be seen on the reflective mirror, you can clean it together with the sensor lens. To do so, use a cotton wool bud slightly (!) moistened in a residue-free solution of water and alcohol to rub carefully over the reflective mirror and, if required, the sensor lens. To clean the housing, wipe it over with a cloth.

19

Do not immerse the device in water or expose it to direct moisture. Avoid the reflective mirror steaming up.

We recommend that you retain the packaging of the NUK Flash Thermometer so that the device can be packed properly again if it needs transporting.

When the device is not in use, we recommend switching it off via OFF mode (see section 3.3.2 OFF mode). Store the device in a cool, dry place out of direct sunlight.

20

Only operate and store the device in accordance with the specified ambient conditions.

3.3 DISPOSAL

The packaging of the NUK Flash Thermometer should be disposed of in an environmentally-sound manner using your local recycling systems.

21

Do not dispose of used batteries in normal waste. Take them to an in-store battery collection system or to your nearest recycling centre for special refuse.


22

Do not dispose of the device as normal waste. Take it to an official collection centre for electrical devices.

4 TROUBLESHOOTING

Message	Problem	Possible cause	Solution
	Device not working properly.	Software error	Remove the battery. Wait one minute and then insert the battery back in (see section 3.1). If the error message appears again, contact your retailer.
	It is not possible to start a new measurement in forehead mode.  symbol is flashing.	The device is not yet ready for a new measurement.	Wait until the symbol has stopped flashing and goes out. When it has stopped flashing, start your new measurement.
	Room temperature above 40.0 °C (104.0 °F).	The device has not been stored in line with the specified storage conditions.	Put the thermometer in a room with a temperature of between 10.0 °C (50.0 °F) and 40.0 °C (104.0 °F) for at least 30 minutes.
	Room temperature below 10.0 °C (50.0 °F).		
	(1) <i>In forehead mode:</i> measurement temperature above +42.2 °C (108.0 °F) (2) <i>Surface measurement:</i> measurement temperature above +80.0 °C (176.0 °F)	The temperature of the object to be measured is outside the specified measuring range (too hot).	Take a measurement in the specified temperature range. If the device is still malfunctioning, contact your retailer.

Message	Problem	Possible cause	Solution
Lo	(1) In forehead mode: measurement temperature below +34.0 °C (93.2 °F) (2) Surface measurement: measurement temperature below -22.0 °C (-7.6 °F)	The temperature of the object to be measured is outside the specified measuring range (too cold).	Take a measurement in the specified temperature range. If the device is still malfunctioning, contact your retailer.
	The device cannot be switched on and the display is empty.	Insulation strips have not been removed from the battery	Remove the insulation strips from the battery.
		Battery empty	Replace the battery (see section. 3.1).
		Battery inserted incorrectly when changed.	Insert the battery correctly (see section. 3.1).

5 TECHNICAL DATA

Temperature measuring range

Forehead mode:

$\pm 0.2^{\circ}\text{C}$ (0.4°F) within $35.0\sim 42.0^{\circ}\text{C}$
($95.0\sim 107.6^{\circ}\text{F}$),

$\pm 0.3^{\circ}\text{C}$ (0.5°F) for other range
 $\pm 0.3^{\circ}\text{C}$ (0.5°F) within $22.0\sim 42.2^{\circ}\text{C}$
($71.6\sim 108^{\circ}\text{F}$),

Surface temperature mode:

others $\pm 4\%$ or $\pm 2^{\circ}\text{C}$ (4.0F) whichever is greater.

Operating temperature:	+10.0~40.0 °C (50.0~104.0 °F), 15%~85% RH
Storage temperature:	-20.0~+50.0 °C (-4.0~122.0 °F) relative humidity \leq 85 %
Transport temperature:	<70.0 °C (158.0 °F) relative humidity \leq 95 %
Accuracy:	
Forehead mode:	± 0.2 °C (0.4 °F) within 35~42 °C (95.~107.6 °F), ± 0.3 °C (0.5 °F) for other range
Surface temperature mode:	± 0.3 °C (0.5 °F) in range 22.0~42.2 °C (71.6~108.0 °F), otherwise ± 4 % (however at least ± 2.0 °C / 4.0 °F)
Air pressure:	800~1013 hPa
Protection class:	BF
Degree of protection:	IP22
Battery:	1 CR2032 lithium battery
Battery life:	approx. 3,000 measurements (~1 year)
Comply with:	ASTM E1965-98, EN ISO 80601-2-56, IEC/EN60601-1-2(EMC), IEC/EN60601-1(Safety) standards, ISO10993, RoHS.
Clinical Bias:	-1.4~ -1.7 °C
Limits of Agreement:	0.98
Repeatability:	0.20°C
Expected service life:	3 years from initial use

Note

The NUK Flash Thermometer converts forehead temperature into its "oral equivalent" (relative to the results of the clinical trial).

Note

The LOT and serial number of the device are located on the inside of the battery compartment lid.

Note

The NUK Flash Thermometer is calibrated during manufacture. If at any time you doubt the accuracy of temperature measurements, contact your retailer immediately.

Guidance and manufacturer's declaration – electromagnetic emissions

The THixyz series is intended for use in the electromagnetic environment specified below. The customer or the user of the THixyz series should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The THixyz series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The THixyz series is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

The THxyz series is intended for use in the electromagnetic environment specified below. The customer or the user of the THxyz series should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the THxyz series, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

Guidance and manufacturer's declaration – electromagnetic immunity

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the THixyz series is used exceeds the applicable RF compliance level above, the THixyz series should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the THixyz series.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity

The THixyz series is intended for use in the electromagnetic environment specified below. The customer or the user of the THixyz series should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the THixyz series requires continued operation during power mains interruptions, it is recommended that the THixyz series be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

**Recommended separation distances between
portable and mobile RF communications equipment and the
ME EQUIPMENT or ME SYSTEM**

The THixyz series is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the THixyz series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the THixyz series as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

6 WARRANTY

In addition to the seller's warranty obligations arising from the contract to purchase, we as the manufacturers provide a warranty of 24 months from purchase of the device as long as the device is used properly and the instructions for use are followed. Evidence of the date of purchase and the type of device should be provided by means of a receipt. The purchaser's statutory rights are unaffected by this manufacturer warranty.

During the warranty period, we undertake to remedy any defects attributable to poor materials or workmanship. Wear parts are excluded from the warranty. Minor deviations from specified quality which do not affect the value or fitness for use of the device are not covered by the warranty. We are likewise unable to accept the warranty if the defects on the device are due to transport damage for which we are not responsible, or, if they are attributable to misuse or lack of care or if the device has been interfered with by persons not authorized by us to do so.

The warranty service will involve repair, replacement of parts or replacement of the device at our discretion. The performance of services under warranty neither extends nor restarts the warranty period. The warranty period for any spare parts fitted will expire with the warranty period for the device as a whole. Extended or different claims, especially those relating to the making good of damage sustained outside the device, are excluded unless liability is compulsory in law.

We will not accept the costs or risks of transport. If a device is sent in without evidence of the purchase date, it will be treated as a repair case. The device will only be repaired following discussion with the customer. Please keep the article number and the lot number for future reference.

These warranty conditions apply to devices bought in Germany. The warranty conditions applicable in the country in question apply to devices not bought in Germany.

23

No modification of this equipment is allowed.

Service Centres

Addresses for service can be found in "Contact addresses" on page 97.

Article no.: 10.256.380

7 EU DECLARATION OF CONFORMITY

The device meets all applicable European directives, as well as the relevant national laws implementing them. These are detailed in the EU declaration of conformity which can be requested from the manufacturer. The declaration of conformity can be found at www.nuk.de.

The NUK Flash Thermometer meets the requirements of ASTM E1965-98 "Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature" and the standards IEC/EN 60601-1 "Medical Electrical equipment - Part 1: General regulations for basic safety and essential performance" and collateral standard IEC/EN 60601-1-2 "Electromagnetic compatibility. Requirements and tests".

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NUK is a registered trademark of
Radiant Innovation Inc.
Medical Technology Promedt Consulting GmbH
MAPA GmbH
MAPA GmbH, Germany