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## User Manual

Please read the instructions before installing and commissioning. Ensure that you retain the instructions for later usage. If you hand over the device for usage to other people, hand over the user manual with it.

## Table of Contents

1. Precaution .....	33
2. Preparation .....	35
Package Contents .....	35
Part Names & Functions.....	36
3. Key Features.....	37
4. How to use the Microscope .....	38
Turn on/off the Microscope .....	38
Display Icons and Meanings.....	39
Focus & Magnifications .....	40
Observe Specimen Slice .....	42
Capture Photo, Video and Playback .....	42
Delete Files .....	43
Use the Polarizer .....	44
Download to Computer .....	45
Preview on Computer .....	45
TV Output .....	46
Charge Battery.....	47
5. Menu Introductions.....	48
Set Resolution .....	48
Capture Mode.....	49
Loop Recording .....	49
Date and Time .....	50
Date Stamp.....	50
TV Out.....	51
Auto Power Off .....	51
Language.....	52
Format SD .....	52
Default Setting .....	53
6. Specifications .....	54
7. Troubleshooting .....	55

8. PC Software Introductions .....	56
Install the Software .....	56
Start Microscope .....	57
Software functions.....	58
Capture Photo .....	58
Capture Video.....	58
Timer.....	59
Calibration and Measurement .....	60
Calibration .....	60
Measurement.....	64
Change Calibration.....	65
Draw and Text.....	68
Software Interface Instructions.....	69
Main Interface Icons .....	69
Main Menu.....	70
Options.....	71
Capture photo/video.....	72
Measurement Interface Instructions.....	73
9. Care and Maintenance .....	74
10. Declaration of Conformity.....	74
11. Disposal.....	74
12. Contact.....	75

## **Welcome**

Thank you for purchasing the 5-inch Portable Stand Alone LCD Digital Microscope with polarizer.

Please take a moment to read through this manual. Its contents will help you to get the best use of the microscope.

## **1. Precaution**

Before using the LCD Digital Microscope, please ensure that you read and understand the safety precautions described below. Always ensure that the LCD Digital Microscope is operated correctly.

- Never operate the LCD Digital Microscope where operation of such devices is restricted. Improper use creates the risk of serious accident.
- Do not try to open the case of the LCD Digital Microscope or attempt to modify the LCD Digital Microscope in any way. Maintenance and repairs should be performed only by authorised service providers.
- Keep the LCD Digital Microscope away from water and other liquids. Do not use the LCD Digital Microscope with wet hands and never in the rain or snow. Moisture creates the danger of fire and electric shock.

- LCD Digital Microscope may become slightly warm during use. This does not indicate malfunction.
- If the exterior of the LCD Digital Microscope needs cleaning, wipe it with a soft, dry cloth.

### **FCC Notice**


This equipment has been tested and complies with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

## 2. Preparation


### ■ Package Contents

- Digital microscope UltraZoom PRO
- 5-V-USB-power supply
- Mini-USB-cable
- Video-cable, K
- CD-ROM with PC-software
- Calibration chart
- Object holding clamps
- User Manual

## ■ Parts Names & Functions

No.	Functions	No.	Functions
1	Micro-SD card slot	8	Height adjustment
2	Video-out	9	Capture
3	HDMI-out	10	5-inch IPS LCD display
4	Focus tube	11	Dimmer for LEDs in the lens tube
5	Polarizer wheel	12	Digital zoom
6	Charging indicator	13	Dimmer for LEDs in the platform
7	LED light for slice	14	USB/DC-in

- 15. Power on/off
- 16. Setup/Settings
- 17. Up
- 18. Confirm
- 19. Down
- 20. Mode (Photo/Video/Playback)


### 3. Key Features

- Magnification: 10x-300x, up to 1200x by digital zoom
- 5-inch IPS panel, 960x540
- Photo resolution: 12M/10M/8M/5M/3M/2M
- Video resolution: 1080p/720p
- Storage: by Micro-SD card (not included)
- Powered by: built-in rechargeable battery, or power adaptor
- Measurement: by software on computer
- Bundle software: PortableCapture Plus for Windows7/8/10, Mac 10.14 or later. Software available on the included CD and [www.dnt.de](http://www.dnt.de)

## 4. How to use the Microscope

### 4.1 Turn on/off the Microscope

Press on/off button to turn the microscope on or off.

## 4.2 Display Icons & Meanings

Icon	Meaning
	Photo capture mode
	Video capture mode
	Playback mode
001618	Quantity of photos that can be captured
12M	Current resolution
	Micro-SD card
	Battery capacity

**Note:** The icon displays will be off when there is no key pressed for 10s.


### 4.3 Focus & Magnifications

You can get different magnifications by rotating the height adjustment wheel and focus tube. The lower the lens, the higher the magnification will be.

- **Low magnification**

Rotate the height adjustment wheel to move the lens upward, and then rotate the focus tube toward left to increase the resolution.

This magnification rate is suitable for getting an overall observation on objects like coins, stamps, jewelry, electronics, etc.


- **Medium magnification**

Rotate the height adjustment wheel to move the lens downward to stop at a height (see left photo), and then rotate the focus tube toward right to get a fine focus.

You can stop at any height and then rotate the focus tube for a better resolution.

- **High magnification**

Rotate the height adjustment wheel to move the lens downward to stop at the lowest height (see right photo), and then rotate the focus tube toward the right to get a fine focus.

This magnification enables you to see the details of an object.



#### 4.4 Observe Specimen Slice

At a High magnification position, turn on the LED at the bottom platform and put a slice on it to see the details of a specimen slice.


#### 4.5 Capture Photo, Video & Playback

Press the capture button to take a photo or record video. Press the mode button to switch among photo, video and playback.


## 4.6 Delete files

- 1) At play back mode, press  to enter the sub-menu, then confirm delete.

- 2) Protect


At play back mode, press  to enter the sub-menu, then confirm protect pictures from delete.

## 4.7 Use the Polarizer


Polarizer, by removing light reflection, is especially useful for observing objects that have a shiny surface like metal, diamond, jewelry, etc.

By rotating the polarizer wheel by 90-degrees, you can activate or deactivate the polarizing effect.

## 4.8 Download to Computer

- 1) Connect the microscope to a computer with a USB cable.
- 2) Choose Mass Storage to view stored photos and videos.


## 4.9 Preview on Computer

- 1) Take out Micro-SD card from the microscope.
- 2) Connect the microscope to a computer with a USB cable.
- 3) Select PC Camera mode
- 4) Start the PortableCapture Plus software. Please refer to Part II on how to use the microscope with the PortableCapture Plus software.

### Note

When working on a computer, the button functions are no longer active, but most of the functions are available in the PortableCapture Plus software.


## 4.10 Output to External Monitor

1. Insert the 1/8-inch plug on the end of the supplied video cable into the microscope's Video port. An HDMI connection is also possible.
2. Insert the RCA plug on the video cable into one of your television's yellow video input connectors. Refer to your television's manual for more information on locating the connectors or HDMI.
3. Turn on the television. If necessary, switch your television to the appropriate input channel. Items that you magnify appear on the television screen.

## 4.11 Charge Battery

Plug the small USB end of the power adaptor into the USB/DC-in port located on the bottom side of the unit. Then plug the other end of the power adaptor into a standard AC outlet. The battery indicator LED at the side of the unit will be in red to show that the battery is charging.


When the battery is finished charging, the battery indicator LED will be in off. Unplug the unit. The battery status can be found on the display.

## 5. Screen Menu Introductions

You can press  to enter the camera setup, and select the options seen below.

### 5.1 Photo Resolution

At photo mode, press , you can choose to select a resolution for photo taken among 12M/10M/8M/5M/3M.

### 5.2 Video Resolution

At video mode, press , you can choose to select a resolution for video recording between 1080P/720P.

### 5.3 Capture Mode

At photo mode, press , you can choose Capture Mode to setup interval capture. For example, configure the microscope to take a photo every 5 seconds.

### 5.4 Loop Recording

At video mode, Press  to set up the length of each video clip. Off means not to turn on loop recording.

## 5.5 Date/Time

Press  and scroll down to set up system date and time.

## 5.6 Date Stamp

At photo mode, press  to select whether to put a date and time mark when capturing photos and videos.

## 5.7 TV-out Mode

To choose between NTSC and PAL for the TV output of the microscope, press  .

## 5.8 Auto Power Off

Press  and scroll down to choose a time for the microscope to power off automatically without key operation.

## 5.9 Language

Press  to set up OSD language.

## 5.10 Format SD

Press  to format the inserted Micro-SD card.

## 5.11 Default Setting

Press  to reset the microscope configuration to default setting.

## 6. Specifications

**Magnification Ratio:** 10x to 300x, up to 1200x digital

**LCD Display:** 5-inch IPS, 960x540

**Photo resolution:** 12M/10M/8M/5M/3M/2M

**Video resolution:** 1080p/720p

**Focus Range:** Manual focus from 5mm to 70mm

**Polarizer:** Switchable polarizer

**Storage:** Micro-SD card

**Output to external monitor:** HDMI or TV

**Powered by:** Built-in rechargeable battery and power adaptor

**Light source:** LEDs with adjustable brightness

**OSD language:** English, German, Spanish, French, Russian, Japanese, Chinese, etc.

**Bundle software:** PortableCapture Plus for Win/Mac  
Software available on the included CD and on [www.cameradownload.net](http://www.cameradownload.net)

**Product Size:** 165(L)x158(W)x267(H)mm

**Weight:** 660g

## 7. Troubleshooting

Problems	Solution
The unit does not turn on.	Charge the battery.
The screen is blank.	Ensure that the power is on. Charge the battery. Verify that the item you want to magnify is positioned correctly under the lens.
Microscope is connected to a television, but no image appears on the television screen.	Ensure that the television is on. Switch to the appropriate video input channel on your television. Ensure that the video/HDMI cable is securely connected to the microscope video-out port and plugged into the correct video input
The battery indicator LED is not showing red, even though the power adaptor is connected and plugged into a wall socket. -or- The battery is not charging.	on your television. Ensure that the charger is plugged into the Microscope's DC-in port.  Ensure that the outlet has power and is not malfunctioning. If the outlet is connected to a wall switch, make sure that the switch is on.

## 8. PC Software Introduction

### 1. Install the Software on PC (Win/Mac)

Insert the CD into the CD-ROM Drive and then open the CD.

Double click to install **PortableCapture Plus** and the installation wizard will guide you through the process.

#### **Note**

The software **PortableCapture Plus** is available on the included CD and [www.cameradownload.net](http://www.cameradownload.net).

## 2. Start Microscope

Connect your Microscope to your PC USB port, start the

software by clicking the  icon generated on the desktop after installation. You will see the following screen pop up.

When the microscope is disconnected from your PC USB port, the following message will appear.

**No Device detected, please connect your microscope directly to your PC USB port.**

### 3. Software Functions

#### ■ Capture Photo


- Click 'Photo' shortcut icon in the main menu bar.
- Click 'Capture > Photo'.

#### ■ Capture Video


- Click 'Video' shortcut icon in the main menu bar.
- Click 'Capture > Video'.

## ■ Timer (Photo/Video)

Click 'Timer' shortcut icon  in the main menu bar.

\* Under Timer Photo, you can set the start time, interval time, and total picture number.


\* Under Timer Video, you can set the start time, interval time, and total video number.

## ■ Calibration and Measurement

Please complete calibration before measurement.

### Calibration

Display Crossing from 'Options>Crossing'.


Microscope default at crossing 'on'

It's used to check the dial on the calibration ruler.


1. Focus microscope on the attached calibration ruler


and adjust the magnification till the picture is clear. Take a picture for the selected area when the dial is parallel to the crossing as shown below.


2. Double-click on the picture to enter the measurement interface as seen below.

3. Click Calibration icon , then move the mouse to the picture. Then click on the beginning and ending points across a known value. Now a dialogue will pop up where the exact value of the measurement sample into 'Actual dimension' needs to be entered.


Then the software will automatically calibrate the magnification rate. Click 'OK' and the adjusted magnification rate will show in the upper right corner.

\* Remark: The measured line should be paralleled.

4. Now you can use 'Any angle Line' under the measurement

icon  to measure the value again. If the value you measure by software is nearly the same as the exact value of the dial, it means that the calibration is completed.

## Measurement

Double-click on the picture to enter the measurement interface. Under the 'Measurement' icon , you can choose any desired option to measure.

Available options are:

### 1. Any Angle Line Measurement

Click from one point and drag to the next point, and then release.

### 2. Continuous Line Measurement

The Continuous Line measurement allows you to measure any continuous desired length.

Click and drag from one point to other desired points, then release.

### 3. Radius Circle Measurement

Click and extend to the desired radius, the radius, circumference and area of the circle will show up accordingly.

### 4. Diameter Circle Measurement

Click and extend to the desired diameter, the diameter, circumference and area of the circle will show up accordingly.

### 5. Three Points Angle Measurement


Define any three points angle and click, the angle value

## Change Calibration

If you need to change the magnification rate and adjust focus to capture the picture, calibrate the magnification rate again.


1. Double click to open the picture, click '  > Reset magnification > Reset picture magnification' as shown below. And then start calibration again (for calibration details, please refer to p. 60).

2. Change the magnification rate and capture several pictures, click on the picture to enter the measurement interface. Then click ' > Set magnification > Set picture magnification', now a dialogue will show up; enter the magnification rate into the dialogue and click 'OK'.


interface. Then click ' > Set magnification > Set picture magnification', now a dialogue will show up; enter the magnification rate into the dialogue and click 'OK'.

3. If you want to use the same magnification rate to


capture pictures, click ' > Set magnification > Set capture magnification'. Now all pictures will use the same magnification rate and you can start the measurement process directly.


## Draw and Text

You can draw or add a text on the picture and save for later analysis.

1. Click the icon  to choose any kinds of drawing under the pull-down arrow as shown below.

2. Click the icon  and set the font for the text in the popup dialogue as below.

## Software Interface Instructions

### Main Interface Icons

Capture photo – It's used to capture photo one by one.

Start/Stop record video.

Timer for photo/video – Set the start time, interval time, and total photo/video number.


Full screen mode

Press ESC or double-click to exit.

## Main Menu


### File

**Photos Directory:** set the directory to store photos;

**Videos Directory:** set the directory to store video;

## Options

**Resolution:** change preview size

**Date/Time:** set to display system date/time on preview window

**Language:** choose different OSD language

**Crossing:** display cross on preview window

**Full-screen mode:** enter into full screen. Click ESC or double click to exit.

## Capture


Photo: take a photo

Video: record a video clip

## Measurement Interface Instructions

Save as

Redo

Copy to  
clipboard

Draw

Copy file

Text note

Previous

Measurement

Next

Calibration

Undo

## 9. Care and Maintenance

Only clean the device with a soft, dry linen cloth. In the case of heavy soiling, it can be slightly damp. Do not clean the device with cleaning agents containing solvents. The device has to be carefully dried with a cloth. Make sure that no humidity gets inside the device. Do not press the display when cleaning!

## 10. Declaration of Conformity

dnt Innovation GmbH, Maiburger Straße 29, 26789 Leer, Germany, hereby declares that the device

„dnt digital microscope UltraZoom PRO“

is in accordance with the essential requirements and other relevant provisions of Directive 2014/30/EU. The declaration of conformity can be found at the following link: [www.dnt.de](http://www.dnt.de)

## 11. Disposal

Do not dispose this device through domestic waste.

To dispose of electronic devices, please refer to the directive of waste for electrical and electronic equipment via the local collection points for old electronic devices.

## 12. Contact

Do you have any questions about the product? Our technical customer service will be happy to provide you comprehensive and qualified information:

Email: [support@dnt.de](mailto:support@dnt.de)

---

1. Edition English 07/2020

Documentation © 2020 dnt Innovation GmbH

All rights reserved. No part of this operating manual may be reproduced in any form without the written permission of the publisher. The operating instructions may have some typographical defects or misprints. However, the information in this operating manual is being checked regularly, and corrections will be done for the next edition. We don't have liability for technical or printing errors and their consequences. All trademarks and copyrights are acknowledged. Changes due to technical progress can be made without notice.

DNT000006-07/2020, Version1.0