TOOLCRAFT

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GB Table of Contents

| | | |
|-----|--|----|
| 1 | Introduction | 23 |
| 2 | Operating Instructions for download | 23 |
| 3 | Intended use | 23 |
| 4 | Delivery contents | 24 |
| 5 | Description of symbols | 24 |
| 6 | Safety instructions | 24 |
| 6.1 | General | 25 |
| 6.2 | Handling | 25 |
| 6.3 | Operating environment | 25 |
| 6.4 | Operation | 25 |
| 6.5 | Power adapter | 26 |
| 6.6 | Li-ion battery | 27 |
| 7 | Product overview | 28 |
| 7.1 | Product | 28 |
| 7.2 | Display elements & control panel buttons | 29 |
| 7.3 | Navigating the menu | 29 |
| 8 | Getting started | 30 |
| 8.1 | Charge the rechargeable battery | 30 |
| 8.2 | Insert the microSD card | 30 |
| 8.3 | Format the MicroSD card | 30 |
| 8.4 | Set the date and time | 31 |
| 9 | Operation | 31 |
| 9.1 | Focus and zoom | 31 |
| 9.2 | Polarizer | 31 |
| 9.3 | Image capture | 32 |
| 9.4 | Video record | 32 |
| 9.5 | Image / video playback | 32 |
| 9.6 | Image and video (delete / protect / slideshow) | 33 |

| | | |
|------|--|----|
| 9.7 | Transfer files to a computer..... | 33 |
| 10 | Settings menu..... | 33 |
| 11 | Computer operation..... | 34 |
| 12 | Cleaning and care..... | 34 |
| 13 | Disposal..... | 35 |
| 13.1 | Product..... | 35 |
| 13.2 | (Rechargeable) batteries..... | 36 |
| 14 | Technical data..... | 36 |
| 14.1 | Product..... | 36 |
| 14.2 | Power adaptor..... | 37 |
| 14.3 | Data according to EU regulation 2019/1782..... | 38 |

1 Introduction

Thank you for purchasing this product.

If there are any technical questions, please contact:

www.conrad.com/contact

2 Operating Instructions for download

Use the link www.conrad.com/downloads (alternatively scan the QR code) to download the complete operating instructions (or new/current versions if available). Follow the instructions on the web page.

3 Intended use

The product is a digital microscope with polarizer that can capture video and still images of viewed objects. Use the product to observe various objects such as coins, stamps, rocks, insects, plants, textiles, circuit boards, etc.

The product is intended for indoor use only. Do not use it outdoors.

Contact with moisture must be avoided under all circumstances.

If you use the product for purposes other than those described, the product may be damaged.

Improper use can result in short circuits, fires, or other hazards.

The product complies with the statutory national and European requirements.

For safety and approval purposes, you must not rebuild and/or modify the product.

Read the operating instructions carefully and store them in a safe place. Make this product available to third parties only together with the operating instructions.

All company names and product names are trademarks of their respective owners. All rights reserved.

The terms HDMI, HDMI High-Definition Multimedia Interface, HDMI Trade Dress and the HDMI Logos are trademarks or registered trademarks of HDMI Licensing Administrator, Inc.

USB4®, USB Type-C® and USB-C® are registered trademarks of USB Implementers Forum.

4 Delivery contents

- Digital microscope
- Power adaptor
- USB-A to USB-C® cable
- HDMI cable (male to male)
- Calibration ruler
- Operating instructions

5 Description of symbols

The following symbols are on the product/appliance or are used in the text:

The symbol warns of hazards that can lead to personal injury.

The symbol warns of dangerous voltage that can lead to personal injury by electric shock.

This product must only be used in dry, enclosed indoor areas. It must not become damp or wet.

6 Safety instructions

Read the operating instructions carefully and especially observe the safety information. If you do not follow the safety instructions and information on proper handling, we assume no liability for any resulting personal injury or damage to property. Such cases will invalidate the warranty/guarantee.

6.1 General

- The product is not a toy. Keep it out of the reach of children and pets.
- Do not leave packaging material lying around carelessly. This may become dangerous playing material for children.
- If you have questions which remain unanswered by this information product, contact our technical support service or other technical personnel.
- Maintenance, modifications and repairs must only be completed by a technician or an authorised repair centre.

6.2 Handling

- Handle the product carefully. Jolts, impacts or a fall even from a low height can damage the product.

6.3 Operating environment

- Do not place the product under any mechanical stress.
- Protect the appliance from extreme temperatures, strong jolts, flammable gases, steam and solvents.
- Protect the product from high humidity and moisture.
- Protect the product from direct sunlight.

6.4 Operation

- Consult an expert when in doubt about the operation, safety or connection of the product.
- If it is no longer possible to operate the product safely, take it out of operation and protect it from any accidental use. **DO NOT** attempt to repair the product yourself. Safe operation can no longer be guaranteed if the product:
 - is visibly damaged,
 - is no longer working properly,
 - has been stored for extended periods in poor ambient conditions or
 - has been subjected to any serious transport-related stresses.

6.5 Power adapter

Do not modify or repair mains supply components including mains plugs, mains cables, and power supplies. Do not use damaged components. Risk of death by electric shock!


- Connect the appliance to a wall socket that can be accessed easily.
- As power supply, only use the supplied mains adaptor.
- Only connect the power adaptor to a normal mains socket connected to the public supply. Before plugging in the power adaptor, check whether the voltage stated on the power adaptor complies with the voltage of your electricity supplier.
- Never connect or disconnect power adaptors if your hands are wet.
- Never unplug the power adaptor from the mains socket by pulling on the cable; always use the grips on the plug.
- For safety reasons, disconnect the power adaptor from the mains socket during storms.
- Do not touch the power adapter if there are any signs of damage, as this may cause a fatal electric shock! Take the following steps:
 - Switch off the mains voltage to the socket containing the power adapter (switch off the corresponding circuit breaker or remove the safety fuse, and then switch off the corresponding RCD protective switch).
 - Unplug the power adapter from the mains socket.
 - Use a new power adapter of the same design. Do not use the damaged adapter again.
- Ensure that cables are not pinched, kinked or damaged by sharp edges.
- Always lay cables so that nobody can trip over or become entangled in them. This poses a risk of injury.

6.6 Li-ion battery

- Never damage the rechargeable battery. Damaging the casing of the rechargeable battery might cause an explosion or a fire!
- Never short-circuit the contacts of the rechargeable battery. Do not throw the battery or the product into fire. There is a danger of fire and explosion!
- Charge the rechargeable battery regularly, even if you do are not using the product. Due to the rechargeable battery technology being used, you do not need to discharge the rechargeable battery first.
- Never charge the rechargeable battery of the product unattended.
- When charging, place the product on a surface that is not heat-sensitive. It is normal that a certain amount of heat is generated during charging.

7 Product overview

7.1 Product

- | | | | |
|----|---|----|---|
| 1 | LCD display | 2 |  Illumination |
| 3 |  Digital zoom | 4 |  Focus dial |
| 5 | Polarizer | 6 |  Record button |
| 7 |  Transillumination | 8 | 5 V/USB USB-C® power input |
| 9 | MicroSD card slot | 10 | HDMI-OUT HDMI port |
| 11 | Control panel | 12 | Height adjustment |
| 13 | Charge indicator | 14 | Battery compartment |

7.2 Display elements & control panel buttons

| Display / button | Description |
|---|---|
|  | Display: battery status |
|  | Display: microSD card inserted |
|  | Display: video mode |
|  | Display: Photo mode (single frame capture) |
|  | Display: Playback mode |
| 00012/01618 | Display: Photos taken / photos remaining |
| 12 M | Display: photo resolution |
|  | Button: power |
|  | Button: settings |
|  | Button: navigate up |
| OK | Button: select |
|  | Button: navigate down |
|  | Button: select operating mode (single image / video / playback) |

7.3 Navigating the menu

| Button | Description |
|---|-------------------------|
|  | Enter and exit the menu |
|  | Navigate UP / DOWN |
| OK | Confirm selection |

8 Getting started

8.1 Charge the rechargeable battery

Important:

Fully charge the rechargeable battery before using for the first time.

1. Connect USB cable to the power adaptor.
2. Connect the USB cable to the **5V/USB** input on the product.
3. Connect the power adaptor to a mains power outlet.
 - The charging indicator will light up.
 - The indicator will switch off when the battery is fully charged.

8.2 Insert the microSD card

Important:

Do not use force when inserting a microSD card into the slot.

Gently insert the microSD card into the slot until it "clicks" into place.

For microSD card specifications, see section: Technical data.

8.3 Format the MicroSD card

To record images and video, the microSD card must be formatted.

Important:

Formatting erases the microSD card. Back up any important files first.

Preconditions:

- ✓ The microSD card is inserted.
1. Press  to switch the power ON.
 2. Press  to enter the menu.
 3. Use  /  to navigate to Format.
 4. Press **OK** and follow the onscreen instructions to format the microSD card.
 5. Press  to exit the menu.


8.4 Set the date and time

1. Press  to switch the power ON.
2. Press  to enter the menu.
3. Press  /  to navigate to Date/Time.
4. Press **OK** to enter settings.
5. Press  /  select values. Press **M** confirm entries.
6. Press  to exit the menu.

9 Operation

9.1 Focus and zoom

1. Place the object to be examined on the platform.
2. Adjust the distance (see: **1**) between the object and the lens.
3. Bring the object into sharp focus (see: **2**).
4. Zoom in further using the digital zoom **-** / **+**.

9.2 Polarizer

- Use the polarizer to remove unwanted light reflections from an object being observed.
- Rotate the polarizer bezel as far as it will go.

9.3 Image capture

Preconditions:

- ✓ A formatted microSD card is inserted.
- 1. Bring the object into focus.
- 2. Press **[M]** until  shows on the display.
- 3. On the product, press  to capture an image.
- The image is saved to the microSD card.

9.4 Video record

Preconditions:

- ✓ A formatted microSD card is inserted.
- 1. Bring the object into focus.
- 2. Press **[M]** until  shows on the display.
- 3. On the product, press  to start recording a video.
 - The record symbol  flashes on the display.
- 4. Press  to stop recording.
- The video is saved to the microSD card.

9.5 Image / video playback

Preconditions:

- ✓ MicroSD card is inserted.
- 1. Press **[M]** until playback mode  shows on the display.
- 2. Use **▲ / ▼** to select a recording.
 - Press **OK** to play or pause videos.
 - Press **▲ / ▼** to adjust the playback speed.
- 3. Press  to enter the playback menu for further options such as: Delete Current, Delete All.

9.6 Image and video (delete / protect / slideshow)

Preconditions:

- ✓ MicroSD card is inserted.
- 1. Press **M** until playback mode  shows on the display.
- 2. Use **▲** / **▼** select an image or video.
- 3. Press  to enter the playback menu and make a selection.

| Selection | Description |
|----------------|--|
| Delete Current | Delete current image or delete all images. |
| Delete All | |

9.7 Transfer files to a computer

Connect the microscope to a computer to transfer captured images and videos from the microSD card to the computer.

Preconditions:

- ✓ A microSD card is inserted.
- 1. Connect the USB cable to the microscope and the other end to a USB port on the computer
- 2. Switch the microscope ON.
 - The microSD card will show up as a drive on the computer.
- 3. Transfer images and videos to the computer.
- 4. After the transfer is complete, eject the drive through the operating system and then remove the USB cable from the computer.

10 Settings menu

| Setting | Description |
|------------------|--|
| Photo Resolution | Set the photo resolution. |
| Time Lapse Photo | Set an interval for image capture (e.g., every 5 seconds). |

| Setting | Description |
|-----------------|---|
| Record Cycle | <ul style="list-style-type: none"> ■ Limit the length of each video clip (e.g., 3 minutes). ■ Recording will stop when the memory card is full. |
| Time Stamp | Enable / disable time stamp on photos. |
| Reticle Switch | Enable / disable reticle overlay (crosshair) on image. |
| Date/Time | Set date, time, and date format. |
| Language | Select language. |
| Format | Format the microSD card. |
| Auto Power Off | Set auto-off time. |
| Default Setting | Reset all settings to the factory default settings. |
| About | View the firmware version. |

11 Computer operation

The product can be operated in combination with a computer and computer software. Download the software and software instructions from:

www.conrad.com/downloads .

12 Cleaning and care

Important:

- Do not use aggressive cleaning agents, rubbing alcohol or other chemical solutions. They damage the housing and can cause the product to malfunction.
- Do not immerse the product in water.
- Avoid touching the lens with any substance.

1. Disconnect the product from the power supply.
2. Clean the product with a dry, fibre-free cloth.

13 Disposal

13.1 Product

This symbol must appear on any electrical and electronic equipment placed on the EU market. This symbol indicates that this device should not be disposed of as unsorted municipal waste at the end of its service life.

Owners of WEEE (Waste from Electrical and Electronic Equipment) shall dispose of it separately from unsorted municipal waste. Spent batteries and accumulators, which are not enclosed by the WEEE, as well as lamps that can be removed from the WEEE in a non-destructive manner, must be removed by end users from the WEEE in a non-destructive manner before it is handed over to a collection point.

Distributors of electrical and electronic equipment are legally obliged to provide free take-back of waste. Conrad provides the following return options **free of charge** (more details on our website):

- in our Conrad offices
- at the Conrad collection points
- at the collection points of public waste management authorities or the collection points set up by manufacturers or distributors within the meaning of the ElektroG

End users are responsible for deleting personal data from the WEEE to be disposed of.

It should be noted that different obligations about the return or recycling of WEEE may apply in countries outside of Germany.

13.2 (Rechargeable) batteries

Remove batteries/rechargeable batteries, if any, and dispose of them separately from the product. According to the Battery Directive, end users are legally obliged to return all spent batteries/rechargeable batteries; they must not be disposed of in the normal household waste.

Batteries/rechargeable batteries containing hazardous substances are labelled with this symbol to indicate that disposal in household waste is forbidden. The abbreviations for heavy metals in batteries are:

Cd = Cadmium, Hg = Mercury, Pb = Lead (name on (rechargeable) batteries, e.g. below the trash icon on the left).

Used (rechargeable) batteries can be returned to collection points in your municipality, our stores or wherever (rechargeable) batteries are sold. You thus fulfil your statutory obligations and contribute to environmental protection.

Batteries/rechargeable batteries that are disposed of should be protected against short circuit and their exposed terminals should be covered completely with insulating tape before disposal. Even empty batteries/rechargeable batteries can contain residual energy that may cause them to swell, burst, catch fire or explode in the event of a short circuit.

14 Technical data

14.1 Product

| | |
|----------------------------|---|
| Input voltage | 5 V/DC 2 A via USB-C® |
| Rechargeable battery | 3.7 V 9.25 Wh, 2500 mAh Li-ion |
| Charging time | 4.5 hrs |
| Operation time | 2 hrs |
| Focus range..... | 5 -70 mm |
| Polarizer | yes |
| Time lapse | 0:00 - 600:60 (mins:secs) |
| Display..... | LCD 178 mm (7") IPS, 1024 x 600 |
| Image capture..... | 16 MP, 12 MP, 10 MP, 8 MP, 5 MP, 3 MP, 2 MP |

| | |
|------------------------------|---|
| Video record | 1440 P 2560 x 1440 @ 30 FPS |
| Video output | HDMI-OUT: 1080i @ 30 FPS |
| Magnification | Optical: 20x to 300x Digital: 300x to 1500x |
| Storage | MicroSD class 10, max. 128 GB File format: FAT32 |
| OSD languages | English, Chinese simplified, Polish, Swedish, German, Italian, Spanish, Slovak, Czech, Por- tuguese, Finnish, Danish, Norwegian, Russian, French, Japanese, Korean |
| Supported OS | Windows® 7/8/10; MacOS 10.6 (or above) |
| Operating temperature | 0 to +40 °C |
| Operating humidity | 20 - 80 % RH (non-condensing) |
| Storage temperature | -5 to +40 °C |
| Storage humidity | 20 - 80 % RH (non-condensing) |
| Dimensions (W x H x D) | 180 x 280 x 165 mm |
| Weight (approx.) | 930 g |

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All rights including translation reserved. Reproduction by any method (e.g. photocopying, microfilming or the capture in electronic data processing systems) requires prior written approval from the editor. Reprinting, also in part, is prohibited. This publication reflects the technical status at the time of printing.

Copyright by Conrad Electronic SE

*3135689_V3_0626_dh_mh_en 18014399982997899-2 I3/O3 en