
MT06MD
User Manual

- English
- Español
- Deutsch
- Francais
- Русский
- 日本語
- 한국의
- Română
- 简体中文

Warranty Details

All NITECORE® products are warranted for quality. Any defective / malfunctioning NITECORE® product can be repaired free of charge for a period of 60 months (5 years) from the date of purchase. Beyond 60 months (5 years), a limited warranty applies, covering the cost of labor and maintenance, but not the cost of accessories or replacement parts.

The warranty is nullified in all of the following situations:

1. The product(s) is/are broken down, reconstructed and/or modified by unauthorized parties.
2. The product(s) is/are damaged through improper use.
3. The product(s) is/are damaged by leakage of batteries.

For the latest information on NITECORE® products and services, please contact a local NITECORE® distributor or send an email to service@nitecore.com

※ All images, text and statements specified herein this user manual are for reference purpose only. Should any discrepancy occurs between this manual and information specified on www.nitecore.com. Sysmax Industry Co., Ltd. reserves the rights to interpret and amend the content of this document at any time without prior notice.

SYSMAX Industry Co. Ltd

TEL: +86-20-83862000
E-mail: info@nitecore.com
Address: Rm1401-03, Glorious Tower, 850 East Dongfeng Road, Guangzhou, China 510600

FAX: +86-20-83882723
Web: www.nitecore.com

(English) MT06MD
User Manual

Features

- Designed for use as a medical inspection light
- Utilises a Nichia 219B LED (CRI>90), max output of 180 lumens
- Photobiologically safe (RG0), colour temperature CCT=5000K
- Turns on in "low" by default, convenience and ideal for pupil inspections
- Powered by 2xAAA batteries, max runtime up to 35 hours
- 3 brightness levels available
- Reverse clicky tail switch provides easy one-handed operation
- High efficiency constant current circuit provides stable and consistent output
- Stainless steel clip
- Aerospace grade aluminium alloy construction
- HAIII hard anodising
- IPX8 (2 meters submersible)
- Suitable for alcohol sterilisation
- Impact resistant to 1.5 meters

Dimensions

Length: 5.04" (128mm)
Head diameter: 0.55" (14mm)
Tail diameter: 0.55" (14mm)
Weight: 0.78oz (22g) (without battery)

Accessories

Spare o-ring

Battery Options

	TYPE	Nominal voltage	Compatible
Primary AAA battery	AAA	1.5V	Y (Recommended)
Rechargeable AAA battery	AAA	1.2V	Y (Recommended)
LiFePO4 battery	10400	3.2V	N (Banned)
Rechargeable Li-Ion battery	10400	3.7V	N (Banned)

Brightness & Runtime


FL1 STANDARD	HIGH	MID	LOW
	180 LUMENS	45 LUMENS	4 LUMENS
	45min	5h15min	35h
	58m (Beam Distance)		
	860cd (Peak Beam Intensity)		
	1.5m (Impact Resistant)		
	IPX8,2m (Waterproof AND Submersible)		

NOTICE: The above data has been measured in accordance with the international flashlight testing standards ANSI/NEMA FL1 using 2x1.2V 750mAh AAA batteries under laboratory conditions. The data may vary slightly during real-world use due to battery type, individual usage habits and environmental factors.

Operating Instructions

Battery Installation

Please install 2xAAA batteries as illustrated.

WARNING: Insert the battery with the positive pole (+) pointing towards the head. The MT06 will not operate with incorrectly inserted batteries

On/Off

With the light turned off, press the switch until a click is heard to enter low output directly. Press the switch again to turn the light off.

Brightness Adjustment

With the light turned on, turning the MT06MD off and on within 2 seconds or simply press the switch halfway to down cycle brightness through "low-mid-high".

The MT06MD will turn on in "low" by default after it is turned off for 2 seconds.

NOTE: The MT06MD will automatically step down in lumens after 5 minutes of use to prevent overheating and extend battery life.

Changing Batteries

Batteries should be replaced when the output appears to be dim or the flashlight becomes unresponsive.

Maintenance

Every 6 months, wipe the threads with a clean cloth followed by a thin coating of silicon-based lubricant.

Warranty Details

All NITECORE® products are warranted for quality. Any defective / malfunctioning NITECORE® product can be repaired free of charge for a period of 60 months (5 years) from the date of purchase. Beyond 60 months (5 years), a limited warranty applies, covering the cost of labor and maintenance, but not the cost of accessories or replacement parts.

The warranty is nullified in all of the following situations:

1. The product(s) is/are broken down, reconstructed and/or modified by unauthorized parties.
2. The product(s) is/are damaged through improper use.
3. The product(s) is/are damaged by leakage of batteries.

For the latest information on NITECORE® products and services, please contact a local NITECORE® distributor or send an email to service@nitecore.com

※ All images, text and statements specified herein this user manual are for reference purpose only. Should any discrepancy occurs between this manual and information specified on www.nitecore.com. Sysmax Industry Co., Ltd. reserves the rights to interpret and amend the content of this document at any time without prior notice.

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