ACE Q

Atemalkoholtester

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Breathalyzer

Operating instructions – EN

Alcotest

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Alcohol Tester

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Test alcolico

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Alcoholímetro

Manual de instrucciones – ES

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Brugsanvisning – DK

Andningsapparat

Bruksanvisning – SE

Alkomat

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Nefesten alkol testleri

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Atemalkoholtester ACE Q

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Operating instructions – EN

Author's liability

The information contained herein has been compiled by ACE Handels- und Entwicklungs GmbH (ACE) conscientiously and on the assumption that it is correct. However, the compiler cannot accept any liability for any use of the information and any consequences arising therefrom. The author is not liable for any costs or expenses arising from the use of the information contained in this document. This information is subject to change without notice.

Introduction

Thank you for purchasing our ACE Q breathalyzer! This device is an easy-to-use breathalyzer with an electrochemical sensor. The device measures the alcohol concentration in the breath and calculates the blood alcohol concentration. It may only be used in accordance with the instructions in the user manual.

This device provides very reliable test results with a measurement tolerance of $\pm 0.05\text{‰}$ (per mille = g/L) at a reference value of 0.50 per mille and 25°C room temperature.

The measuring unit ‰ (per mille) = basic setting in the EU.

The measurement results determined in per mille are displayed with 3 decimal places (e.g.: 0.21%) and can be interpreted 1:1 (e.g. the display shows the value 0.36 = 0.36 per mille). Measured values in % BAC are displayed with 4 digits.

The electrochemical sensor integrated in this breathalyzer must be calibrated. To maintain the measuring accuracy of this device and the sensor warranty, we recommend regular calibration. You can send the device for calibration to the address found on www.ace-technik.com for calibration.

The device may only be used in accordance with the instructions in the operating instructions. All readings obtained from

this breathalyzer are for guidance only and are not legally binding. Breath tests should be carried out using clip-on mouthpieces (a new mouthpiece should be used for each test for reasons of hygiene and to avoid falsified readings).

Note on country-specific measuring units

If necessary, please inform yourself about the current country-specific alcohol limits, which are available, for example, on the website of the Ministry of Transport or from automobile clubs.

Contents

- 1 ACE Q digital breathalyzer with electrochemical sensor
- 6 mouthpieces
- 2 AAA alkaline batteries (1.5 V)
- 1 storage bag
- 1 instruction manual (multilingual)
- 1 calibration return slip

Designation of the individual components

To conduct a test with the most accurate result possible:


- A) Wait at least 15 minutes after the last sip of alcohol until the alcohol test. Do not eat, drink or smoke during this time.
- B) Do not open the housing.
- C) Do not store the breathalyzer in a cold, hot, damp or poorly ventilated place.
- D) Keep the breathalyzer out of the reach of children.
- E) After a longer storage period, carry out three "run-in tests" before carrying out the actual alcohol test.
- F) Do not drop the breathalyzer.
- G) Replace the battery if its capacity is low. Reduced battery capacity leads to incorrect test results.
- H) Do not carry out a test in noisy (sound waves interfere with the sensor), draughty or poorly ventilated areas.
- I) After a test result of 1.00 per mille or higher, the sensor requires a longer recovery time. In this case, wait at least three minutes until the next test.
- J) Do not carry out a test in a very cold, very hot or unventilated, smoky environment.
- K) Store the breathalyzer in such a way as to prevent the ingress of dust and liquids.
- L) Do not blow too hard into the device, as this will damage the sensor (especially if a lot of saliva gets into the device when blowing).

Preparation for the test

1. Install 2 AAA batteries (alkaline, 1.5 V) in the battery compartment. To do this, carefully remove the battery compartment cover on the back of the device and insert the two batteries. Note the battery polarity (+/- pole).
2. Insert a mouthpiece into the side opening for attaching the mouthpieces. Make sure that the mouthpiece is firmly attached. Use a new mouthpiece for each test subject (for reasons of hygiene and to maintain measurement accuracy, as a used mouthpiece may contain residual alcohol inside).

Instructions for conducting a test

1. Press the power button for 2 seconds to switch the device on. The preset alarm limit as well as the date and time appear.

- DATE / TIME
25-07-30 09:31
- NEXT CALIBRATION
363 DAYS

2. The appliance warms up, after which the words "START BLOW" appear on the display.
3. Blow constantly for approx. 5 seconds without stopping until you hear two beeps. Tip: It is best to blow as if you were whistling without a sound.
4. At the end of the analysis, the alcohol concentration, measured in per mille, is displayed for approx. 4 seconds in lock mode; the power button and menu button are deactivated for this period. The screen with the displayed measurement result then remains and you have the following operating options due to the suspension of the button operating lock:

#82 25/04/16 12:30
0.000 ‰
⏻-Retest

Retest: Press the power button briefly to initiate another test.

5. The device switches off automatically after 3 minutes without performing any further tests. You can also switch the device off immediately by simply holding the power button for 3 seconds.

ACE's Recommendations for Proper Testing

- Always use a new mouthpiece for each person testing! ACE is not responsible for any damage or injury that may result from the use of one mouthpiece by several people and for several tests.
- Always use the hygienically and individually packaged ACE mouthpieces. The test device comes into contact with human excretions (sweat, saliva, blood, etc.) both inside and outside during each test.
- Use gloves when removing the mouthpieces used for a test from the device, especially if the test has been carried out by someone else. ACE Technik (www.ace-technik.com) has suitable gloves in its range.

Important Notes

- If you have eaten something, drunk alcohol or smoked a cigarette, wait at least 15 minutes before starting the measurement.
- For measured values below 0.070‰ (OR 0.035 mg/L OR 0.007 % BAC - depending on the setting), the test result is displayed as 0.00 (‰ OR mg/L OR % BAC - depending on the setting). This is because at a low alcohol level, other substances may be present in the breath that influence the result of the breath sample.
- This breathalyzer is a high-quality and sensitive measuring device. In order to maintain the measuring accuracy of the electrochemical sensor, the device should be used regularly. If the device is not used for a longer period of time, the sensor may dry out. We recommend using the breathalyzer at least twice a month for breath tests (even if you have not drunk alcohol).
- The processes described may change due to modifications to the product. This may occur without prior notice and without any adjustment to the operating instructions and does not constitute grounds for complaint, return or price reduction. It is the responsibility of the user to recognize deviations between the operating instructions and the actual operation of this breathalyzer and to react accordingly.

Only the use of original ACE mouthpieces

- To ensure the accuracy of the measurement results.
- To guarantee 100% health safety.
- To be eligible for the full device warranty

WARNING

Even if the device indicates an alcohol level in the "safe" range, the indication is not a license to drive a vehicle or perform any other activity that should not be performed under the influence of alcohol.

Explanation of the Symbols and Icons

- 1. Total number of tests performed
- 2. Current temperature
- 3. Current time
- 4. Remaining battery power

Explanation of error / warning messages

Display	Meaning
BLOW ERROR	If you see "Flow" on the display, the breath sample was not strong enough for the measurement. The device warms up again automatically. If "blow" is displayed, blow hard and long into the mouthpiece again until you hear a signal tone.
LOW BATTERY	If you see this symbol in the display, the battery is low. Please replace the batteries with 2 new AAA batteries.
CALIBRATION DUE DATE HAS PASSED	Calibration is required because the device has exceeded its recommended service life. The device must be calibrated to yield accurate results. Send the device in for calibration.
OVER NUMBER OF TEST ALLOWED	Calibration is required because the device has exceeded the recommended number of tests. The device must be calibrated to yield accurate results. Send the device in for calibration.

BREATH SAMPLE WITH EXCESSIVE SALIVA	Excess saliva entered the sensor during the test. Wait 20 minutes and then try again. Contact your seller's customer service if this error code is displayed for a longer period of time.
OUT OF TEMP RANGE	The test cannot be performed because the appliance is switched on outside the operating temperature. Repeat the test after the appliance has reached the temperature within the operating range.
SAMPLING ERROR	The test cannot be performed due to an issue in the sampling system. Contact your local seller or customer service center for proper repair service.

Calibration

Calibration is not mandatory, but the accuracy of the integrated sensor decreases steadily without calibration. Send the device to the address stated on the calibration return bill. **We recommend sending it as a parcel so that proof of delivery or an inquiry can be made if the device cannot be found.**

The processes described may change due to modifications to the product. This can occur without prior notice and without adapting the operating instructions and is not a reason for complaint, return or price reduction. It is the responsibility of the user to recognize deviations between the operating instructions and the actual operation of the ACE Q and to react accordingly

Notes on the batteries

The device referred to in these operating instructions is supplied with 2 AAA alkaline batteries. These batteries are not intended to be recharged. The use of a rechargeable battery in the device is also prohibited, as this will damage the device. If the battery is empty, it must not be disposed of with household waste. Used batteries may contain harmful substances that can damage the environment and health. Please dispose of batteries/rechargeable batteries at retail outlets or at local authority recycling centers. The return is free of charge and required by law. Please only dispose of discharged batteries in the containers provided and tape the terminals of lithium batteries. Battery

recycling is the easiest environmental protection exercise in the world - please support it!

The waste icons shown below mean that old batteries and rechargeable batteries must not be thrown away with household waste. They can be recycled and valuable raw materials such as iron, zinc and nickel can be recovered.

The abbreviations under the waste garbage cans stand for:

- Hg: Battery contains mercury
- Cd: Battery contains cadmium
- Pb: Battery contains lead

WEEE-REG-Nr. from ACE: 89961185

LUCID reg. no. from ACE: DE4358374179837

Technical data*

Sensor	Electrochemical sensor
Accuracy	+/- 0.05 ‰ at 1.00 ‰ and 25°C (per mille = basic setting) OR +/- 0.005 %BAC at 0.10 % BAC and 25°C OR +/- 0.025 mg/L at 0.50mg/L and 25°C
Indication of the BAC (blood alcohol content) concentration)	0.00 up to 5.00 ‰ (per mille = basic setting) OR 0.000 up to 0.500 % BAC OR 0.00 up to 2.50 mg/L
Warm-up time	5 seconds at 0.50 per mille (the higher the measured per mille value, the longer the sensor needs to evaporate the alcohol content measured in the breath)
Device dimensions	115 x 53 x 25mm (H x W x D)
Weight	90g (including 2 AAA batteries)
Energy supply	2 1.5V AAA alkaline batteries (= energy for up to 1,000 tests)

Temperature conditions	Operating temperature: 5°C to 50°C Storage temperature: 0°C to 50°C
Calibration recommendation	After 12 months or 1,000 tests (whichever is earlier). You can send the device to the service address given at www.ace-technik.com for calibration or purchase a replacement sensor there.

*Specifications are subject to change without notice for functional improvements.

MENU OPTIONS

Accessing the menu:

Press the menu button for two seconds during the warm-up process

Controls in the menu:

- Upwards: Press the power button briefly.
- Downwards: Press the menu button briefly.
- Confirm the selection: Press the power button for 2 seconds.
- Back to the warm-up process: Press the menu button for 2 seconds.

To save the settings:

Press and hold the power and menu button together for 2 seconds until you hear a beep.

TEST RESULT:


The test results can be checked under the "TEST RESULT" tab. The test number, test result, date and time are also saved. The last test result is displayed first.

- To display the previous measurement result, quickly press the power button.
- To display the next measurement result, quickly press the menu button.

- To return to the menu, press the menu button for 2 seconds

ALARM LIMIT

The alarm limit can be set in this menu. By default, this is set to 0.5 ‰ (per mille) or 0.050% BAC (BAC) = same measured value.

The following limits can be set:

- 0.1‰ - 1.0‰ (0.010 ~ 0.100 %BAC)
- Switching off the limit
- To increase the alarm limit, briefly press the Power Button.
- To decrease the alarm limit, briefly press the Menu Button.
- To save and return to the menu list, press and hold the Power Button for 1 second.
- To return to the menu list without saving, press and hold the Menu Button for 1 second.

UNIT

The measuring unit of the device can be set in this menu item.

The options for measuring units are:

- ug/L
- mg/L
- ‰
- mg/100mL
- %BAC
- g/210L
- To change the units, briefly press the Power Button or the Menu Button.
- To save and return to the menu list, press and hold the Power Button for 1 second.
- To return to the menu list without saving, press and hold the Menu Button for 1 second.

BRIGHTNESS

The brightness of the screen can be adjusted in this menu item.


The adjustable value range:

- 5 levels (1=low, 5=high)
- To brighten the display, briefly press the Power Button.
- To darken the display, briefly press the Menu Button.
- To save and return to the menu list, press and hold the Power Button for 1 second.
- To return to the menu list without saving, press and hold the Menu Button for 1 second.

SOUND

The volume of the sound can be adjusted here.


The sound can be set in this menu.

The adjustable value range:

- 5 levels (1=low, 5=high)
- To increase the volume, briefly press the Power Button.
- To decrease the volume, briefly press the Menu Button.
- To save and return to the menu list, press and hold the Power Button for 1 second.
- To return to the menu list, press and hold the Power Button for 1 second.
- To return to the menu list without saving, press and hold the Menu Button for 1 second.

INFORMATION

Further information on the device can be found here.

The following information is listed in the menu item:

- Firmware version
- Last calibration date
- The total number of all tests

LANGUAGE

Set the language

The settable options include:

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- To change the language, briefly press the Power Button or the Menu Button.
- To save and return to the menu list, press and hold the Power Button for 1 second.
- To return to the menu list without saving, press and hold the Menu Button for 1 second.

PIN-CODE

The correct PIN code must be entered in this menu to access the advanced menu option.

The default PIN code is 0000, and the PIN code can be set from 0000 to 9999 in the advanced menu option.

- To change to higher digits of PIN code, press and hold the Power Button for 1 second.
- To increase the number of PIN code, briefly press the Power Button.
- To decrease the number of PIN code, briefly press the Menu Button.
- To enter advanced menu options, press and hold the Power Button for 1 second at the first digit (thousand) of PIN Code and select YES by pressing and holding the Power Button for 1 second for confirmation.
- To return to the menu list, press and hold the Menu Button for 1 second.

Once you have entered the correct code, you will be taken to the **ADVANCED MENU**:

ADDITIONAL MENU OPTIONS:

TEST LIMIT

The maximum permitted number of measurements can be set in this menu.

The adjustable value range:

- 50 ~ 2000 tests
- OFF

- To increase the test limit, briefly press the Power Button.
- To decrease the test limit, briefly press the Menu Button.
- To save and return to the advanced menu list, press and hold the Power Button for 1 second.
- To return to the advanced menu list without saving, press and hold the Menu Button for 1 second.

SET TIME


The device time and date can be set in this menu.

- Device time is displayed as yyyy/mm/dd, hh:mm:ss.
- To change digits of time and date, press and hold the Power Button for 1 second.
- To increase the time and date, briefly press the Power Button.
- To decrease the time and date, briefly press the Menu Button.
- To save the mode change, press and hold the Menu Button for 1 second or press the Power Button at the last digit of time (seconds). Select YES by pressing and holding the Power Button for 1 second for confirmation.
- If no changes were made, you can return to the menu list, by pressing and holding the Menu Button for 1 second.

DATE LIMIT

The date for the next calibration can be set in this menu.

The adjustable value range: 1 ~ 24 months

- To increase the date limit, briefly press the Power Button.
- To decrease the date limit, briefly press the Menu Button.
- To save and return to the advanced menu list, press and hold the Power Button for 1 second.
- To return to the advanced menu list without saving, press the Menu Button for 1 second.

CHANGE PIN

The PIN code can be set in this menu.

The adjustable value range:

- 0000 ~ 9999
- To change to higher digits of PIN code, press and hold the Power Button for 1 second.
- To increase the number of PIN code, briefly press the Power Button.
- To decrease the number of PIN code, briefly press the Menu Button.
- To save and return to the advanced menu list, press and hold the Power Button for 1 second at the first digit (thousand) of PIN code or press and hold the Menu Button for 1 second at any digit. Select YES by pressing and holding the Power Button for 1 second for confirmation.
- To return to the advanced menu list without saving when no changes have been applied, press the Menu Button for 1 second.

Warranty and liability

ACE Handels- & Entwicklungs GmbH (ACE) guarantees the product to which this guarantee and liability declaration is attached - in all versions and designations - for two years after purchase against defects in material and workmanship; the sensor is guaranteed for one year. This warranty includes the electronics and the housing or refers to the sensor without sampling components. The warranty does not cover components that are inevitably destroyed during operation (e.g. fuses). The warranty is valid from the date of delivery or from the date of first commissioning (whichever is shorter), if this is specified in the literature accompanying the ACE product.

WARRANTY AND LIMITATION OF LIABILITY

1. Product restriction

ACE disclaims all other warranties, express or implied, including the implied warranties of merchantability and fitness for a particular purpose. In particular, the warranty does not apply to calibration. Calibration is mandatory after sensor replacement. If a product requires sensor replacement and/or calibration, the user may request a cost estimate from ACE. If the user rejects the quote, ACE will return the product unrepaired/uncalibrated, and the product will be labeled accordingly as such. If the product does not fulfill this warranty, the purchaser's claim and ACE's obligations are limited exclusively and at ACE's discretion to the replacement or repair of the affected goods or to a refund of the purchase price for the rejected goods.

2. Restriction of use

In no event shall ACE be liable for any additional, incidental or consequential damages, including, without limitation, lost profits and loss of use, arising out of the sale, manufacture or use of the products sold under this disclaimer, whether based on contract, negligence or other tortious action.

The product to which this guarantee and liability declaration is attached, irrespective of the version and designation, cannot be used in a legally binding manner to determine the fitness of a person to drive a vehicle or to carry out an action that could endanger themselves or others under the influence of alcohol.

Mandatory requirements for full product warranty include: careful inspection of all products on receipt for damage, calibration suitable for the purchaser's particular application and proper use, and repair and maintenance in accordance with the instructions in the product literature supplied by ACE. Repair or maintenance work carried out by unqualified personnel will invalidate the warranty, as will the use of non-approved consumables or spare parts. As with all technically demanding products, it is a fundamental and mandatory requirement of ACE's warranty that all users of these products are familiar with the use, capabilities and limitations of the product as described in the relevant product literature.

3. Restriction on consultation

The buyer confirms that he/she has verified the intended use and suitability of the purchased goods for any purpose under its sole responsibility. It is expressly agreed between the parties that any technical or other advice given by ACE regarding the use of the goods or services has been given free of charge and at the buyer's sole risk. Therefore, ACE shall have no obligation or liability whatsoever arising from the recommendations or their results.

DEVICE CARE TIPS

Important notes on calibration

Any breathalyzer purchased from ACE can also be calibrated by ACE. Of course, we also calibrate devices purchased elsewhere, provided these devices are identical in construction to the devices sold by ACE.

You can see the due date of the next calibration on the calibration sticker on your breathalyzer. If the sticker is missing, then the calibration should take place 12 months after the purchase date: e.g. purchase date: 13.05. = next calibration on 13.05. and then on 13.05. of the following year.

You can send the device to the service address given at www.ace-technik.com for calibration.

Why is calibration necessary?

Calibration is not mandatory, but the accuracy of the integrated sensor will steadily decrease without calibration and the warranty will also expire if the calibration intervals are not observed.

Factors influencing the measurement accuracy:

- The more often the breathalyzer is used, the more the zero point shifts.
- At high alcohol concentrations, the zero point shifts more strongly.
- The sensitivity of the sensor also decreases with devices that are used little or not at all.

ACE recommends testing with the breathalyzer at least once (1x) a week, regardless of whether you have drunk alcohol. The device should also be calibrated regularly to maintain measurement accuracy.

Please also note that various factors influence the measurement result, including but not limited to:

- Alcohol quantity
- Period of alcohol consumption
- Gender
- Age
- Weight
- Body size
- Drinking habit
- Health condition
- Vitality of the liver
- Size of the lung
- Intake of food during alcohol consumption, and
- Influence of medication.

Important notes on carrying out the test and self-assessment

You should wait at least 15 minutes after the last sip of alcohol, as there may still be residual alcohol in the oral cavity, which affects the test result. Smoking should also be avoided before a test, as cigarette smoke damages the sensor.

A new mouthpiece should be used for each test - not only for hygienic reasons, but also because a mouthpiece that has already been used can influence the next test result. There could be residual alcohol in the mouthpiece, which could influence the new test result if the same mouthpiece is used again (e.g. this could indicate alcohol in the blood of a sober person).

If "Flow" is shown on the display, the breath sample was not sufficient, or it was ended too early. During your breath test, blow constantly for approx. 4 seconds as if you were whistling without a sound. The breath test is finished as soon as you hear a double beep.

If the same person carries out 2 tests in succession, they should wait at least 3 minutes between the two tests, as the measured values may be too different if the waiting time is too short. The air exhaled from the lungs after a test is no longer enriched with the same alcohol content as when the test was carried out shortly before. Therefore, an immediate second test will show a lower reading.

All measured values determined with this breathalyzer are only a guide value and are not legally binding.

Why the perceived alcohol level and the measurement result often do not (or cannot) match

In an alcohol test, there are many different factors that influence the measured value. In addition to gender, age, height and body weight, these also include the body fat percentage and the amount of fluid in the body

In addition, the general physical and mental state or the respective daily condition of the subject must also be taken into account. The type of meal before the test (type of contents and stomach filling level), and the intake of medication can also influence the measured value.

Other factors that influence alcohol metabolization include the type of drink (percentage of alcohol content in the drink as well as its carbon dioxide and sugar content) and the rate of alcohol consumption.

Some people quickly have an increased alcohol content in their blood, for others it takes longer until a drunken state is noticeable. Around 40-50% of the alcohol consumed is in the body within 15 minutes, after 30 minutes it is 65% and after 70 minutes at the latest, 100% of the alcohol previously consumed has been absorbed by the body. However, it should also be noted that some people, patquikly break down the alcohol content in their blood, whereas others still have an increased alcohol content in their blood even after a longer period of time.

The EU Declaration of Conformity for this product can be downloaded from www.ace-technik.com.

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Bu rn iin AB Uygunluk Beyanı www.ace-technik.com adresinden indirilebilir.