



Content

- General information 3
- Introduction 4
- Display & operating buttons 5
- Preparing before operating 6
- Turn on & measure 6
- Calibration 8
- Changing scale & temperature unit 10
- Turning off 11
- Cleaning & maintenance 11
- Disposal 11
- Technical data 12
- Error codes 12
- Models and scales 13

⚠ Carefully read through the operating manual even if you have prior experience with KERN refractometers.

1. General information

1.1 Intended use

The refractometer is a measuring instrument for determining the refractive index of transparent substances in liquid or in some cases also in the solid state. It is used to observe the behaviour of light as it passes from a prism with known properties to the substance being tested. Use of the refractometer for other purposes is contrary to its intended use and may be hazardous. The manufacturer shall not be liable for any damages caused by improper use.

1.2 Warranty


The warranty shall be void in the event of:
• Failure to observe the instructions in the operating manual
• Use for purposes other than those described
• Modifications or opening the device housing
• Mechanical damage and/or damage resulting from media, liquids, natural wear and tear



This digital refractometer cannot measure any liquid that is highly corrosive to metal or glass. When measuring liquids that are corrosive to plastics or react chemically with plastics, be careful not to drop the measured liquid onto the shell. Otherwise it will corrode the shell.

2. Introduction

2.1 Description




2.2 Scope of delivery

1x Storage box | 1x Digital refractometer | 1x Operating manual | 1x AAA Battery 1.5 V | 1x Pipette | 1x Screwdriver

3. Display & operating buttons

3.1 Description display & operating buttons




2.2 Scope of delivery

1x Storage box | 1x Digital refractometer | 1x Operating manual | 1x AAA Battery 1.5 V | 1x Pipette | 1x Screwdriver


4. Preparing before operating

4.1 Install the battery

Turn the screw counterclockwise to open the battery hatch.




5. Turn on & measure




5. Turn on & measure

5.1 Turn on

The multi-function display shows the current scale number




5.1 Turn on




5. Turn on & measure

5.2 Measure

After turning on, clean the sample tank with distilled water and then dry it. Now fill the sample up to the mark, close the cover and press "READ".




5.2 Measure




5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.



5.3 Average value measurement




5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.




5.3 Average value measurement




5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.




5.3 Average value measurement




5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.




5.3 Average value measurement




5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.




5.3 Average value measurement




5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.



5.3 Average value measurement



5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.



5.3 Average value measurement



If the calibration was not completed successfully, an error code appears in the display. Here, for example, A01.



Further error codes can be found in the appendix.




We recommend calibrating the refractometer,

- when commissioning
- after a strong shock
- after longer transport
- after a change of location with a large temperature difference
- if the device has not been used for a long time


Always use distilled water and make sure that the refractometer, the water and the environment are at the same temperature.

6. Calibration

6.1 Changing scale




6.2 Changing temperature unit




7. Changing scale & temperature unit

7.1 Changing scale



7.2 Changing temperature unit



8. Turning off

If without any operations for 1 minute, the instrument would be automatically turned off.

9. Cleaning & maintenance

1. To avoid damages to the prism and the sample tank, clean them with distilled water after each use.

2. Dry it with a soft cloth afterwards.

3. Do not use hard or abrasive objects for cleaning.

4. Do not leave any residue in the sample tank.

5. If the refractometer is not going to be used for a longer time, remove the battery and store it at a cool and dry place.

10. Disposal

The packaging consists of environmentally friendly materials which can be disposed of via local recycling facilities.

The device and storage box should be disposed

11. Technical data

Scale + accuracy + resolution	Depends to the model
Temperature	0,0 – 40,0 °C / 32,0 – 104,0 °F
Automatic Temperature Compensation	Yes
Minimum sample volume	0.2 - 0.3 ml (Marking ring)
AUTO-OFF	60 seconds
Averaging measurement	15 measurements
Battery	1 x AAA 1.5 V
Lifetime of the battery	Approx. 10.000 measurements
Overall dimensions LxWxH	125x65x30 mm
Net weight	140 g (without battery)

12. Error codes

code	Instructions
A01	Beyond the scope of calibration temperature. (0.0°C~40.0°C)
A02	During calibration, no solution or solution wrong.
A03	This instrument has a hardware failure.

13. Models and scales

Model	Scale	No.	Range	Unit	Resolution	Accuracy
ORM 508M	Refractive Index	502	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1RS	Brix	501	0.0-50.0	%	0.1%	±0.2%
ORM 15U	Refractive Index	502	1.330-1.5177	nD	0.0001nD	±0.0003nD
ORM 25U	Fructose	501	0.0-68.9	%	0.1%	±0.2%
ORM 1H0	Glucose	502	0.0-59.9	%	0.1%	±0.2%
ORM 1H0	Brix	501	0.0-68.9	%	0.1%	±0.2%
ORM 1H0	Refractive Index	502	1.3330-1.5177	nD	0.0001nD	±0.0003nD
ORM 1H0	Lactose	501	0.0-16.5	%	0.1%	±0.2%
ORM 1H0	Maltose	502	0.0-15.6	%	0.1%	±0.2%
ORM 1H0	Dextran	503	0.0-15.6	%	0.1%	±0.2%
ORM 1H0	Brix	504	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Honey Water	501	5.0-38.0	%	0.1%	±0.2%
ORM 1H0	Honey Saturated	502	30.0-48.0	%	0.1%	±0.2%
ORM 1H0	Brix	503	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	504	1.3330-1.5177	nD	0.0001nD	±0.0003nD
ORM 1H0	Salinity NaCl (%)	501	0.0-28.0	%	0.1%	±0.2%
ORM 1H0	Salinity NaCl (‰)	502	0-280	‰	1‰	±2‰
ORM 1H0	Specific Weight	503	1.000-1.220	-	0.001	±0.002
ORM 1H0	Refractive Index	504	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Salinity Seawater	501	0-100	‰	1‰	±2‰
ORM 1H0	Chlorinity Seawater	502	0-57	‰	1‰	±2‰
ORM 1H0	Specific Weight	503	1.000-1.070	-	0.001	±0.002
ORM 1H0	Brix	504	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	505	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Alcohol Malt.	501	0-72	%	1%	±1%
ORM 1H0	Alcohol Vol.	502	0-80	%	1%	±1%
ORM 1H0	Brix	503	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Plo	505	0.0-30.5	°P	0.1	±0.3
ORM 1H0	SG Wort	506	1.000-1.130	-	0.001	±0.002
ORM 1H0	Brix	507	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	508	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Oechsle	509	0.0-150	°Oe	1	±2
ORM 1H0	Vol%	510	0.0-22.0	%	0.1%	±0.2%
ORM 1H0	KMW (Babo)	511	0.0-25.0	°	0.1	±0.2
ORM 1H0	Brix	512	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Oechsle France	513	0-250	°Oe	1	±2
ORM 1H0	Vol%	514	0.0-22.0	%	0.1%	±0.2%
ORM 1H0	KMW (Babo)	515	0.0-25.0	°	0.1	±0.2
ORM 1H0	Brix	516	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	517	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	518	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	519	0.00-25.00	-	0.01	±0.20
ORM 1H0	Brix	520	0.00-20.00	%	0.01%	±0.20%
ORM 1H0	Refractive Index	521	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	522	0.0-12.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	523	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	524	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Urine Cal	525	1.000-1.060	-	0.001	±0.002
ORM 1H0	Urine Dog	526	1.000-1.060	-	0.001	±0.002
ORM 1H0	Brix	527	0.0-50.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	528	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	529	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	530	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	531	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	532	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	533	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	534	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	535	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	536	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	537	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	538	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	539	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	540	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	541	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	542	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	543	0.0-51.0	%	0.1%	±0.2%
ORM 1H0	Refractive Index	544	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1H0	Refractive Index	545	0.0-51.0	%	0.1%	±0.2%