ION USB SCALE

health

Quickstart Guide

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
Guida rapida

Snelstartgids

www.ionhealthonline.com

MANUAL VERSION 1.2

ENGLISH

BOX CONTENTS: USB Scale, USB flash drive, 4 AAA batteries, Quickstart Guide

CAUTION

- Before starting any exercise program, we recommend consulting your physician for a physical exam and to discuss your program.
- The scale is not water-resistant. Do not submerge it in liquid or expose it to wet conditions.
- Avoid rough use or severe impact to the scale.

IMPORTANT!:

- Your bare feet must touch the sensors for your readings to be stored.
- Use the included USB flash drive to store/transfer readings. Other USB flash drives may not be compatible with USB Scale.

FEATURES

1. **SET** – Turns scale on (if it is off). Press this button when the scale is on to select the User. Use it as an "Enter" button to move to the next screen in the LCD.
2. **LEFT/RIGHT** – Decreases/increases the value of the current setting or toggles through available options in the LCD.
3. **SENSORS** – These sensors help calculate body fat percentage, water percentage, muscle mass percentage, and bone mass percentage (based on your body height, age, weight, etc. in your User profile). Your bare feet must touch these sensors for the scale to store these readings.
4. **LCD** – Displays information related to the scale's operation.
5. **USB PORT** – Connect the included USB flash drive to this port. As you weigh yourself, the data will be stored in this drive (your bare feet must touch the sensors for your readings to be stored). When you want to load the data into your software, connect the drive to a free USB port on your computer. (See TRANSFERRING RESULTS TO A USB FLASH DRIVE.)
6. **UNIT (not shown in diagram)** – This button on the underside of the scale is used to set or reset various settings. See INITIAL SETUP for more information.

INITIAL SETUP

1. Open the battery compartment on the back of the scale. Insert 4 AAA batteries and close the compartment. (Be sure the "+" and "-" terminals are facing the correct way.)
2. The number flashing on the LCD represents the year. Use the LEFT/RIGHT buttons to adjust the year and press SET to confirm.
3. Repeat Step 2 for the Month, Day, Hour, and Minute.
4. Press UNIT (on the underside of the scale) once to switch it on. Press it again (or multiple times, if necessary) to toggle through the available measurement units – "kg" (kilograms), "st/lb" (stones/pounds), or "lb" (pounds).

USER PROFILE SETUP

1. Press SET to turn the scale on (if it is off). The user number "P1," "P2," etc. will flash.
2. Use the LEFT/RIGHT buttons to select the desired user (P1-P8). Press SET to confirm.
3. Use the LEFT/RIGHT buttons to select the gender (male or female) and athleticism (athletic male or athletic female, represented by an additional "runner" icon). Press SET to confirm.
Note: Select an "athletic" profile if you are close to your recommended weight for your height and age, as verified with your physician.
4. Use the LEFT/RIGHT buttons to select your height. Press SET to confirm.
5. Use the LEFT/RIGHT buttons to select your age. Press SET to confirm.

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ENGLISH

WEIGHING YOURSELF

- Place the scale on a flat, hard, level surface. The scale may not produce an accurate reading when placed on soft surfaces, like carpets, mats, etc.
- If you are the only person who will normally use the scale, skip to Step #3. If you need to change to a different user profile, press SET so the "P__" is flashing, and use the LEFT/RIGHT buttons to select the correct user.
Note: You can also step on the scale immediately without selecting a user, which will power it on and take a reading. However, it will either (1) automatically save the reading to a user profile with a similar reading (± 0.3 kg / 6.6 lbs. of your current weight reading) or (2) ask you to manually select which user profile to save it to (if there are multiple profiles with readings similar to yours) – wait until the desired profile number ("P__") appears, then press either the LEFT or RIGHT button to select the left or right profile number shown in the display.
- Step on the scale and remain still. When the weight unit measurement (lb, kg, or st/lb) appears, that is your current weight.

- If you are wearing socks, shoes, etc.,** it will measure your weight only. The reading will not be stored.
- If you weigh yourself barefoot (recommended),** several bars will appear on the LCD and gradually deplete, indicating it is taking all your other measurements. Remain still until all of them are gone. The LCD will then show your weight, body fat percentage ("FAT"), water mass percentage ("TBW"), muscle mass percentage ("MUS"), and bone mass percentage ("BONE"), as well as the user profile to which it is saving the reading. The measurements will be shown in sequence three times before the scale turns off. The reading will be stored in USB Scale's internal memory.

If there is no movement and no buttons are pressed, the scale will shut off automatically after 10 seconds.

TRANSFERRING RESULTS TO A USB FLASH DRIVE

- Make sure the scale is off or in a normal mode (i.e., you are not setting your user information).
- Insert a USB flash drive into the USB PORT on the side of the scale. Use the included USB flash drive to store/transfer readings. Other USB flash drives may not be compatible with USB Scale.
Note: *Do not keep the USB flash drive inserted when not in use – it can drain the scale's battery.*
- Press and hold SET for 2-5 seconds. The user number ("P1," "P2," etc.) will flash on the LCD.
- Use the LEFT/RIGHT buttons to select the user whose data you want to transfer.
- Press SET to transfer the data. **"MEM"** will flash on the LCD. The scale will turn off when the transfer is complete.
Note: *Once a user's results are transferred, the scale's memory for that user is erased. (The files are named "CYUSER01," "CYUSER02," etc. and are stored on root directory of the USB flash drive.)*
- Remove the USB flash drive.
- Connect the USB flash drive to a free USB port on your computer. *To use the included software to view and track your data, double-click the "USB Scale Software" folder, double-click "SetupEN.exe," and follow the on-screen instructions to install it. See the software manual PDF to learn how to use the software.*

OTHER MESSAGES

In some cases, the LCD may display messages prompting you to do something:

- MEM** – Indicates that some action is occurring with USB Scale's memory. This is not an error message. (If this icon flashes when transferring readings to a USB flash drive, it means there are no readings left for that user profile in the scale's internal memory.)
- Lo** – Low battery. Replace the batteries.
- Weight overload. Remove the weight on the scale.
- USB icon (flashing)** – Insert a USB flash drive.
- USB icon (solid)** – A USB flash drive is connected.
- FULL (without a USB icon)** – The scale's memory is full. Transfer the data onto a USB flash drive.
- FULL + USB icon (solid)** – Your USB flash drive's memory is full. Transfer the data to your computer.
- USB icon + Err** – The scale may not be able to read or write to the USB flash drive. Try reformatting the drive or use a different one.

SPECIFICATIONS

Weight units: lb, kg, st/lb

Measurements: weight, body fat %, body water % (TBW), muscle mass % (MUS), bone mass %

Power: 4 AAA alkaline batteries

Number of user profiles: 8

Age Range: 10-80

Height Range: 3' 3.5" – 7' 2.5" (100-220 cm)

Max capacity: 397 lbs. / 180 kg / 28 st, 5 lb

Warranty and technical support information can be found online at www.ionhealthonline.com or by calling **401-658-3743** (U.S.) or **01252 896 080** (U.K.).

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