D. Körperanalyse

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B. Preparation

1. Remove battery strip.

Note: Immediately (as long as „0.0“ is still being displayed) place the body analysis scale on a level surface and wait until the scale powers down automatically.

Then begin data input. Otherwise, it might be possible that the weight of the scale would erroneously be added to the initial weighting process.

The same might occur during regular operation if the scale is held in the hand prior to the weighting process.

2. For all measurements, place scale level on a hard surface (not carpet).
3. Cleaning and Care: clean only using slightly damp cloth. Do not use solvents or abrasives. Never submerge scale in water.
4. A country-specific option to switch from kg/cm to lb/in or st/in is possible by pushing the key on the bottom of the scale to switch units.

Attention! Slipping hazard in case of wet surface.

C. Data input

For the correct functionality of the body analysis scale it is necessary to enter the personal data of the respective person and then immediately step on the scale with bare feet.

1. Place scale on a table to enter data.
2. Start data input (push SET key for 1 second)
3. Select memory slot (P1 ... P8) (< = minus, > = plus) and confirm (SET).
4. Enter body height (< = minus, > = plus) and confirm (SET).
5. Enter age (< = minus, > = plus) and confirm (SET).
6. Select gender using </> ( = male,  = female) and confirm (SET).

7. Set the value of your individual activity level according to the table (< = minus, > = plus) and confirm (SET).

Activity level*

1-3	Persons with average activity level, with light work while sitting or standing	Max. 2 hours/day standing or moving
4	Athletes, people with plenty of movement	Intensive training at least 5 hrs/week
5	Heavy manual laborer, highly active athletes	Intensive training at least 10 hrs/week

8. The number of the memory slot is being displayed.

Push SET key and hold for **3 seconds** until "0.0" appears on the screen.

9. Now immediately place the body analysis scale on the floor, wait until the screen displays „0.0“ and then step on the scale with bare feet.

The initial measurement is indispensable and serves to store the data for later automatic recognition of persons. If no person steps on the scale, data input must be repeated.

The first body analysis is performed after a person steps on the scale.

This scale is equipped with a Time-Out mode. It will automatically switch off after approx. 40 seconds if no keys are pressed.

Data must be entered again after battery change.

***Note:** For this model, the selection of the activity levels 1-3 does not result in different measurement results.

D. Body analysis

For most precise results, please step on the scale while unclothed.

1. Step on the scale. The weight and - for the duration of the analysis (5 to 10 seconds) - a moving bar indicator will be displayed. Please stand still until the analysis is complete.
2. After automatic recognition of persons, the scale will display the body weight, body fat content (☉), body water content (≡), and muscle mass (♣) along with the corresponding status bar.

The average range for the respective person is located in the median 3 segments of the bar indicator.

3. The scale will then display the memory slot and one or more of the following symbols. This means that one or more body values are not within average range.

kg Body weight
 ☉ Body fat content
 ≡ Body water content
 ♣ Muscle mass

If the scale symbol (⚠) is displayed, your analysis values are within average range!

4. Step off the scale. After displaying the memory slot, the scale switches off automatically.

In rare cases it is not possible to automatically allocate the analysis results. In this case, the scale will alternately display the memory slots of persons coming closest to the analysis result. Since weight measurement is already completed, you may now step off the scale and use the > or < key (< key for the first memory slot displayed, > key for the second memory slot displayed) to confirm your individual memory slot. Following this, body weight, body fat and body water content, as well as muscle mass will be displayed.

In case of weight differences of more than +/- 3 kg compared to the last measurement, the user will not be recognized, which means that data input must be repeated.

Please note that when stepping on the scale while wearing footwear or socks, only the weight is being determined and a body analysis is not performed or performed incorrectly!

You will achieve continuously precise results if you always stand in the same position on the scale if possible.

E. Messages

1. Batteries spent.
Insert new batteries
2. Overload: Over 150 kg.
3. Implausible analysis result: review personal programming
4. Improper foot contact: clean scale or feet and ensure sufficient moisture on skin. If skin is too dry, slightly wet feet or perform measurement preferably after bath or shower.

F. Technical data

Load capacity x graduation =
 max. 150 kg x 100 g
 Graduation of body fat content: 0.1%
 Graduation of body water content: 0.1%
 Graduation of muscle mass: 0.1%
 Body height: 100-250 cm
 Age: 17-99 years
 8 individual memory slots
 Batteries: 2 x 3 V CR2032

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International

08:30 bis 12:00 Uhr

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For 3 years as of the date of purchase, Soehnle guarantees the removal of defects due to material or processing errors, free of charge, by way of repair or replacement. In case of a guarantee claim, please return your scale to your seller, together with purchase receipt and guarantee card.

**EC Conformity**

This device complies with the applicable EC Directive 2004/108/EC.

Battery disposal EC Directive 2008/12/EC

Pb - Cd - Hg

Batteries are not a part of your regular household waste. You must return batteries to your municipality's public collection or wherever batteries of the respective type are being sold.

Batteries containing hazardous materials are marked with these symbols:

Pb = Battery contains lead

Cd = Battery contains cadmium

Hg = Battery contains mercury

Disposal of electric and electronic devices EC Directive 2002/96/EC

This product is not to be treated as regular household waste but must be returned to a collection point for recycling electric and electronic devices. Further information is available from your municipality, your municipality's waste disposal services, or the retailer where you purchased your product.

Consumer service

If you have questions or comments, the following points of contact are available:

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