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**Facial tanner**

**Bronzeur facial**

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I/B Version  
**120515**

**Typ 834**

230 V~50 Hz 75 W

**Typ 836**

230 V~50 Hz 105 W

**Typ 838**

230 V~50 Hz 150 W

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**efbe-Schott®**

# User Instructions

**Note:** Read and follow the User Instructions carefully before using the tanning system. Keep the User Instructions handy.

## What you should know about skin tanning

### The Sun

is the source of life on earth. Everyone knows from the change of seasons and personal experience during vacations that the light and heat of the sun affect our health and well-being in general. If your skin is bronzed, you look more attractive and healthy and feel better. Skin tan also had a biological function in that it protects your skin from sunburn, e.g. before you go on a sun vacation.

### Tanning

is caused by a certain part of the light spectrum. Tanning does not depend on whether the light is natural sunlight or artificial light.

The radiation of the sun contains a visible portion, which we perceive as light and an invisible portion, which we perceive, for example, as heat. The light also contains radiation that makes the skin turn red or even burns it. Radiation that gives the skin a tan is invisible radiation, viz. ultraviolet (UV) radiation consisting of UVA, UVB and UVC radiation. Intensive UVA radiation with a very low proportion of UVB and no UVC produces the best tan. The earth's atmosphere keeps the UVC radiation away, whereas UVB and UVA radiation reach the earth's surface and tan the skin. However, as the UVB proportion of solar light is rather high, pale skin is irritated easily, becomes red and sun-burned.

### The UVA tanning system

has better proportions of UVA and UVB radiation than natural light. The tanning UVA radiation is much more intense than in sunlight whereas the harmful UVB proportion is only a fraction of that of the solar proportion.

A person using the tanning system usually gets a suntan quicker than from exposure to sunlight and the solarium does not burn the skin. The tan is similar to suntan. However, everybody's skin may respond differently to exposure to sunlight.

### Skin type and exposure time

We all know how our skin responds when it is exposed to direct sunlight. A similar reaction can be expected from exposure to artificial light. Those whose skin can stand much sunlight will bronze quickly; those whose skin is not tanned by the sun will not profit much from the UVA tanning system.

## Skin type by reaction of the unexposed skin to natural sunlight

### Skin Type    Skin reaction and ethnic affiliation

I	Always gets sunburn quickly, hardly any or only little tan, even after repeated exposure (Celtic type)
II	Gets sunburn almost always; moderate tan after repeated exposures (fair-skinned, European type)
III	Gets sunburn moderately often; progressing tanning after repeated exposures (dark-skinned, European type)
IV	Gets sunburn rarely; quick and strong tan (Mediterranean type)

The use of a tanning system is only recommended to persons fitting the description of skin type II, III or IV. Skin-type I persons are unable to develop a normal tan either from the sun or exposure to artificial light. Exposure, if at all, should be after consulting a doctor.

Decide your exposure time after careful and considering your experience with natural light exposure.

The following exposure sequences are recommended for the various skin types.

### Typ 834 - Radiation distance: 5 cm

#### Exposure schedule

Skin type		II	III	IV	
1st exposure	The time between the exposures must be 48 hours minimum	10	10	10	min.
2nd exposure		11	12	13	min.
3rd exposure		13	15	16	min.
4th exposure		14	17	20	min.
5th exposure		16	19	23	min.
6th exposure		17	21	26	min.
7th exposure		18	24	29	min.
8th exposure		20	26	32	min.
9th exposure		21	28	36	min.
10th exposure		22	31	39	min.
11th exposure		24	33	42	min.
12th exposure		25	35	45	min.
First exposure period: Initial dose		10 (= 100 J/m <sup>2</sup> )	10 (= 100 J/m <sup>2</sup> )	10 (= 100 J/m <sup>2</sup> )	min.
Maximum period of exposures: Final dose		25 (= 250 J/m <sup>2</sup> )	35 (= 350 J/m <sup>2</sup> )	45 (= 450 J/m <sup>2</sup> )	min.
Maximum number of exposures per year:		60	43	33	
Yearly dose		15	15	15	kJ/m <sup>2</sup>

### Typ 836 - Radiation distance: 5 cm

#### Exposure schedule

Skin type		II	III	IV	
1st exposure	The time between the exposures must be 48 hours minimum	7	7	7	min.
2nd exposure		8	9	9	min.
3rd exposure		9	10	12	min.
4th exposure		10	12	14	min.
5th exposure		11	13	16	min.
6th exposure		12	15	18	min.
7th exposure		13	17	21	min.
8th exposure		14	18	23	min.
9th exposure		15	20	25	min.
10th exposure		16	21	27	min.
11th exposure		17	23	30	min.
12th exposure		18	25	32	min.


First exposure period: Initial dose	7 (= 100 J/m <sup>2</sup> )	7 (= 100 J/m <sup>2</sup> )	7 (= 100 J/m <sup>2</sup> )	min.
Maximum period of exposures: Final dose	18 (= 250 J/m <sup>2</sup> )	25 (= 350 J/m <sup>2</sup> )	32 (= 450 J/m <sup>2</sup> )	min.
Maximum number of exposures per year:	60	43	33	
Yearly dose	15	15	15	kJ/m <sup>2</sup>

## Typ 838 - Radiation distance: 10 cm

### Exposure schedule

Skin type		II	III	IV	
1st exposure	The time between the exposures must be 48 hours minimum	6	6	6	min.
2nd exposure		7	8	9	min.
3rd exposure		8	9	11	min.
4th exposure		9	11	13	min.
5th exposure		10	12	15	min.
6th exposure		11	14	17	min.
7th exposure		12	15	19	min.
8th exposure		13	17	21	min.
9th exposure		14	18	23	min.
10th exposure		14	20	25	min.
11th exposure		15	21	27	min.
12th exposure		16	23	29	min.
First exposure period: Initial dose		6 (= 100 J/m <sup>2</sup> )	6 (= 100 J/m <sup>2</sup> )	6 (= 100 J/m <sup>2</sup> )	min.
Maximum period of exposures: Final dose		16 (= 250 J/m <sup>2</sup> )	23 (= 350 J/m <sup>2</sup> )	29 (= 450 J/m <sup>2</sup> )	min.
Maximum number of exposures per year:		60	43	33	
Yearly dose		15	15	15	kJ/m <sup>2</sup>

## SUN BOX Typ 834 836 838


## Operating position

Set the unit on a table of suitable height.

## Starting the tanning system

Make the mains connection. Turn the timer knob clockwise to set the exposure time. This starts the UVA tanning system and the timer.

The system stops automatically at the end of the exposure time.


## Memo slide

The tanning effect can be controlled by setting the memo slide.

Every user uses on of 4 memo slides.

Starting at 1, the slide is set to the next number at the end of each exposure.

## Hints for the exposure

For exposures of the face or the upper part of the body choose a seat height that the head is opposite the upper half of the radiation area. The lamps will be warm only a little and could be touched riskless during using the unit. An even tanning of the face will be realized by a little turning your head to the left and right side.

## Important notes:

- Supply only 230-volt alternating current to the appliance.
- Do not cover the vent slots at the back of the appliance while in use.
- For reasons of electric safety, do not use water on the appliance or immerse it in water.
- Do not use the appliance in damp rooms or near a sink.
- Do not use the appliance in an explosion-hazard area.
- Do not use the appliance with a defective timer.
- Only store the appliance after it has cooled down (cooling period ca. 15 minutes).
- This appliance can be used by children over 8 years old and people with reduced physical, sensory or mental capacities or with a lack of experiences or knowledge if they have been given supervision or instructions concerning the safe use of the appliance and the risks involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children unless they are older than 8 years and supervised.
- Keep the appliance and its cord out reach of children less than 8 years, electric appliances. Never let a child use an electric appliance without supervision.
- It is possible that the lights flicker on new devices or show moving shadow effects. This is for the browning function without meaning.
- Certain materials that change colour in sunlight (e.g. paintings) should not be exposed to the UV radiation of the tanning system for a prolonged time.
- Unplug the appliance after use.
- Non-users, especially children, should not be present when the appliance is in use.

- UV appliances should **not** be used by:
    - People under 18;
    - People with the tendency to get freckles;
    - People with natural red hair colour;
    - People with atypically decolourised skin areas;
    - People who have more than 16 moles (2 mm diameter or bigger) on their body;
    - People who have atypical moles (atypical moles are described as asymmetric moles with a diameter that is bigger than 5 mm with different pigmentation and irregular borders);
    - People who suffer from sunburn;
    - People who cannot tan at all or who can not tan without getting sunburn when they are exposed to the sun;
    - People who sunburn easily when they are exposed to the sun;
    - People who have had a history of repeated, heavy sunburns during childhood;
    - People who suffer or suffered from skin cancer or who are predisposed to it;
    - People whose relatives of the first degree have had black skin cancer (melanoma);
  - People who are in medical treatment because of suffering that comes along with photosensitivity;
  - People who get photosensitising medication;
- The exposures should not exceed the personal minimal erythral dose (MED), i.e. the minimal UV dose that causes perceptible reddening of the skin. If an erythema (reddening of the skin) occurs hours after an exposure, no further exposures should take place. After a week the exposures can be taken up again starting from the beginning of the exposure plan.
  - If unexpected effects occur, such as itching within 48 h after the first exposure, ask your doctor for advice before further exposures.
  - People (including children) who are not able to use the appliance in a safe way, due to their physical, sensorial or mental capacity or their lack of experience or knowledge, must never use the appliance without supervision or instructions of a responsible person.

## **Please note the following information**

- Ultraviolet radiation of the sun or from artificial light sources can cause skin or eye damage. The biological effect depends on the type and amount of irradiation and also on the sensitivity of the skin exposed to it.
- Excessive irradiation can cause sunburn. Overly frequent exposure to ultraviolet radiation of the sun or an artificial light source can lead to early ageing of the skin and increased risk of skin tumours. These risks increase according to the accumulated UV radiation. Exposure at a young age increases the risk of skin damage later on in life.
- The surface of the unprotected eye can be inflamed and in certain cases excessive exposure to radiation can cause retinal damage. Cataract can be the consequence of many repeated irradiation exposures.
- In case of pronounced individual sensitivity or allergic reactions to ultraviolet radiation, ask your doctor for advice before the start of the exposure.

## **The following precautionary measures have to be respected:**

- Always use the supplied safety goggles, contact lenses and sunglasses are no surrogate for the safety goggles;
- Remove cosmetics before exposure and do not use any sun protecting agents or self-tanning lotion;
- Certain medical conditions or side effects of medication can be heightened because of ultraviolet radiation. In case of doubt, consult your doctor;
- Let an interval of at least 48 hours between the first two exposures;

- When using the appliance, do not sunbathe the same day;
- Follow the recommended exposure time, exposure intervals and keep the recommended distance to the tubes;
- Consult your doctor if your skin starts swelling, wounds or pigmented spots appear;
- Sensitive skin areas such as scars, tattoos and private parts need to be protected before exposure.
- Consult your doctor if your skin starts swelling, wounds or pigmented spots appear;
- Sensitive skin areas such as scars, tattoos and private parts need to be protected before exposure.

## Cleaning and care

- Pull the mains lead before starting to clean.

## Bulb replacement

The replacement of the lamps is carried out as follows:

- Loosen the fastening bolts at the top of the unit's rear side. Remove the shutter. Turn each lamp a quarter of a turn to remove them from the socket.
- If lamps are suspected defective but work in other sockets, the starter may be the cause. You may also test their function replacing them against each other.

## Caution

UV fluorescent bulbs have to be replaced by bulbs of the same type.

## Changing the starters

Turn the starter counterclockwise at the unit's rear panel. Insert the new starter turning the starter clockwise.

## Limited warranty

We are neither responsible for broken glass of fluorescent tubes nor for handling damages on defective units sent in with insufficient packing. We recommend to keep the original packing for this purpose.

**Caution!** If fluorescent tubes do not light after switching them on, turn the tubes slightly back and forward.

Our warranty does not cover non-defective units which were sent in with fluorescent lamps which did not light, but are non-defective.

**Caution!** If the mains cord of the units gets damaged, it shall be replaced by the manufacturer or authorised service centres or other qualified people to avoid any hazard.

**Caution!** Should you give away the appliance to another user, make sure to hand over the instruction manual together with the appliance.

## Replacement tubes

Typ 834	4 pieces	UV-A lamps	Typ Cleo	15 W	29 cm
Typ 836	6 pieces	UV-A lamps	Typ Cleo	15 W	29 cm
Typ 838	8 pieces	UV-A lamps	Typ Cleo	15 W	29 cm

Like all fluorescent tubes, tanning system tubes contain specific admixtures. These tubes are classed as hazardous waste and must be disposed of according to law.

**If you order spare tubes or starters, please quote the serial data on the tube or starter and the type number of your tanning system.**

Modifications and repairs must be carried out by our after-sales service staff.

Spare parts and customer service available in Germany.

## PROTECTION OF THE ENVIRONMENT- DIRECTIVE 2002/96/EC

In order to preserve our environment and protect human health, the waste electrical and electronic equipment should be disposed of in accordance with specific rules with the implication of both suppliers and users.

For this reason, as indicated by the symbol  on the rating label or on the packaging, your appliance should not be disposed of as unsorted municipal waste. The user has the right to bring it free of charge to a municipal collection point performing waste recovery by means of reuse, recycling or use for other applications in accordance with the directive.

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efbe Elektrogeräte GmbH