efbe-Schott[®]

Infrarotlampe

Infrared lamp

Lampe infrarouge

Infraroodlamp

I/B Version
100906

Art.-Nr. IR 801-818

230V ~

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GB OPERATING INSTRUCTIONS

IMPORTANT SAFEGUARDS

- ❑ Please read these instructions carefully before using the appliance.
- ❑ Only operate the appliance for the intended use as described in the manual.
- ❑ If you hand over the appliance to another user, make sure to join the instruction manual.
- ❑ Check that your mains voltage corresponds to that stated on the appliance.
- ❑ **Danger of burning!** Make sure not to touch the hot surface of the appliance. During use there is rising heat above the appliance. **Caution! Hot surface!**
- ❑ Accessible surfaces (outer housing) can become very hot during use.
- ❑ This appliance can be used by children over 8 years old and people with reduced physical, sensory or mental capacities or with a lack of experiences or knowledge if they have been given supervision or instructions concerning the safe use of the appliance and the risks involved. Children shall not play with the appliance.
- ❑ Cleaning and user maintenance shall not be made by children unless they are older than 8 years and supervised.
- ❑ Keep the appliance and its cord out of reach of children less than 8 years.
- ❑ Never leave the appliance unsupervised when in use. Keep out of reach of children or incompetent persons.
- ❑ From time to time check the cord for damages. Never use the appliance if cord or appliance shows any signs of damage. Should the cord be damaged, it must be replaced by a competent qualified electrician^(*). All repairs should be made by a competent qualified electrician^(*).
- ❑ In order to avoid electric shocks, never immerse the appliance in water or any other liquid for any reason whatsoever. Never place it into the dishwasher.
- ❑ Unplug the appliance before cleaning it.
- ❑ Make sure the appliance has totally cooled down before cleaning or storing it.
- ❑ All repairs should be made by a competent qualified electrician^(*).
- ❑ Never use the appliance outside and always place it in a dry environment.
- ❑ Never move the appliance by pulling the cord. Make sure the cord cannot get caught in any way. Do not wind the cord around the appliance and do not bend it.
- ❑ Stand the appliance on a table or stable surface in order to prevent it from falling down.
- ❑ **Never use the Infrared lamp in the bathroom, the shower or any other humid room.**

^(*) Competent qualified electrician: after sales department of the producer or importer or any person who is qualified, approved and competent to perform this kind of repairs in order to avoid all danger. In any case you should return the appliance to this electrician.

GENERAL INFORMATION

This infrared lamp is effective in the short-wave infrared area and therefore, it is very useful for heat radiation.

The infrared lamp gives full heating immediately after it has been switched on.

To obtain ideal heat radiation, keep a distance of at least 50cm (between the skin and the lamp).

You can perform heat radiation as often as you like. However, insert pauses of several hours between two radiations.

A radiation treatment of approximately 15 minutes is enough in most cases. Avoid too strong falls of temperature after the radiation.

In case you are suffering from serious pain, inflammations or symptoms of unknown nature, please contact your physician before use.

OPERATION

Please put the infrared lamp on a stable non-slip surface before use.

- Remove all packaging materials (carton strips) before using the appliance.
- Plug in the appliance.
- Bring the lamp to the desired radiation position (Caution: remove all inflammable objects and metal objects from the radiation area!).
- Orient the heat radiation towards the area that should be treated. The lamp's orientation can be adjusted vertically in three different positions. When the light is switched on, the temperature around the bulb might become very high. Therefore, when adjusting the orientation of the lamp, only take the lamp by its lateral side.
- Perform the radiation treatment (see general information).
Protect your eyes by closing them during the radiation treatment.
- Unplug the appliance after use and let the appliance cool down for several minutes. Then, insert the cord in the cord storage compartment and store the appliance.

TECHNICAL FEATURES

Typ 801 / 808 / 810 / 818	230 V~ 100 W
Typ 812	230 V~ 150 W

CLEANING

Only clean the appliance with soft damp cloth. Before cleaning, the lamp must be unplugged and cooled down.

SPARE PARTS

Typ 801 / 808 / 810 / 818	Philips R 95 E / E 27 / 100 W
Typ 812	Philips Par 38 E / E 27 / 150 W

PROTECTION OF THE ENVIRONMENT– DIRECTIVE 2002/96/EC

In order to preserve our environment and protect human health, the waste electrical and electronic equipment should be disposed of in accordance with specific rules with the

implication of both suppliers and users. For this reason, as indicated by the symbol  on the rating label or on the packaging, your appliance should not be disposed of as unsorted municipal waste. The user has the right to bring it free of charge to a municipal collection point performing waste recovery by means of reuse, recycling or use for other applications in accordance with the directive.

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efbe Elektrogeräte GmbH