# BRAUN

## ThermoScan®

IRT 3020

<b>GB</b> English .....	4
<b>BG</b> Български .....	14
<b>CZ</b> Česky .....	24
<b>GR</b> Ελληνικά .....	34
<b>HR</b> Hrvatski .....	44
<b>HU</b> Magyar .....	54
<b>PL</b> Polski .....	64
<b>RO</b> Română .....	74
<b>RS</b> Српски .....	84
<b>RU</b> Русский .....	94
<b>SI</b> Slovenščina .....	104
<b>SK</b> Slovenčina .....	114
<b>TR</b> Türkçe .....	124
<b>SA</b> اللغة العربية .....	143

This product is manufactured by Kaz Europe SA under a license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

ThermoScan is a registered trademark of Kaz Europe SA

 Kaz Europe SA  
Place Chauderon 18  
CH-1003 Lausanne, Switzerland

©2011, Kaz Europe SA


www.kaz-europe.com

Designed in Germany Made in PRC

US Patents: 6,022,140 D397,624

IRT3020 EE

PN: 31IM3020205  
25MAY11

## Product Description

- |                   |                         |
|-------------------|-------------------------|
| 1. Protective cap | 5. Lens filter detector |
| 2. Lens filter    | 6. Display              |
| 3. Probe          | 7. Battery door         |
| 4. «start» button |                         |

## English

The Braun ThermoScan thermometer has been carefully developed for accurate, safe and fast temperature measurements in the ear. The shape of the thermometer probe prevents it from being inserted too far into the ear canal which can hurt the eardrum.

**However, as with any thermometer, proper technique is critical to obtaining accurate temperatures.**

**Therefore, read the use instructions carefully and thoroughly.**

### Important safety information

- The operating ambient temperature range for this thermometer is 10 – 40 °C (50 – 104 °F).
- Do not expose the thermometer to temperature extremes (below –20 °C / –4 °F or over 50 °C / 122 °F) or excessive humidity (> 95 % RH).
- This thermometer must only be used with genuine Braun ThermoScan Lens Filters. Other lens filters can lead to inaccuracy.
- To avoid inaccurate measurements always use this thermometer with a new, clean lens filter.
- If the thermometer is accidentally used without a lens filter, clean the lens after use (see «Care and cleaning»).
- Keep lens filters out of reach of children.
- This thermometer is intended for household use only.
- Use of this thermometer is not intended as a substitute for consultation with your physician.
- This thermometer must stay in stable ambient (room) temperature for 30 minutes before operating.

## Body temperature

Normal body temperature is a range. The following table shows that this normal range also varies by site. Therefore, measurements from different sites should not be directly compared.

### Normal temperature ranges of different sites in the body:

Axillary <sup>1</sup> :	34.7 – 37.3 °C	94.5 – 99.1 °F
Oral <sup>1</sup> :	35.5 – 37.5 °C	95.9 – 99.5 °F
Rectal <sup>1</sup> :	36.6 – 38.0 °C	97.9 – 100.4 °F
ThermoScan <sup>2</sup> :	35.8 – 38.0 °C	96.4 – 100.4 °F

- A person's normal temperature range tends to decrease with age. The following table shows normal ThermoScan ranges by age.

### Normal ThermoScan ranges by age<sup>2</sup>:

0 – 2 years	36.4 – 38.0 °C	97.5 – 100.4 °F
3 – 10 years	36.1 – 37.8 °C	97.0 – 100.0 °F
11 – 65 years	35.9 – 37.6 °C	96.6 – 99.7 °F
> 65 years	35.8 – 37.5 °C	96.4 – 99.5 °F

Everyone has their own normal temperature range, and this range also varies throughout the day. It is easy to determine normal temperature range using the Braun ThermoScan. Practice taking temperatures on yourself and healthy family members to determine your normal temperature ranges.

Note: When consulting your physician, communicate that the ThermoScan temperature is a temperature measured in the ear and if possible, quote the individual's normal ThermoScan temperature range as an additional reference.

1. Chamberlain, J.M., Terndrup, T.E., New Light on Ear Thermometer Readings, Contemporary Pediatrics, March 1994.
2. Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, Annals of Emergency Medicine, January 1995.

## How to use your Braun ThermoScan

1. To achieve accurate measurements, make sure a new, clean lens filter (2) is in place before each measurement.

2. Gently squeeze the opposite ends of the thermometer to pull off the probe cap (1). Do not use force in removing the cap.



3. To turn the thermometer on, push the «start» button (4).

During an internal self-check, the display shows all segments.


The thermometer is ready to use after two beeps are sounded and this display appears.

4. Perform an ear tug to straighten the ear canal. This gives the thermometer a clear view of the eardrum.

An ear tug is best performed by using your free hand to grasp the outer edge of the top half of the ear. To take your own temperature, wrap your free hand around the back of your head and tug your ear from behind.

Children under 1 year:  
Pull the ear straight back.

Children aged 1 year to adult:  
Pull the ear up and back.

5. While tugging the ear, fit the probe (3) snugly into the ear canal, then push the «start» button (4).

After one second a beep will signal the end of the measuring process, the result is shown on the display (6).

After 2 short beeps, the ear symbol will stop flashing.

6. For the next measurement, remove the used lens filter and put on a new, clean lens filter by snapping it into position.

Fit the probe snugly into the ear canal, then push the «start» button.

The Braun ThermoScan ear thermometer turns off automatically after 60 seconds of inactivity.

## Temperature taking hints


- Always replace with disposable lens filters to maintain accuracy and hygiene.

- A measurement taken in the right ear may differ from the measurement taken in the left ear. Therefore, always take the temperature in the same ear.
- The ear must be free from obstructions or excess earwax build-up to take an accurate reading.
- External factors may influence ear temperatures, including when an individual has:
  - been lying on one ear or the other
  - had their ears covered
  - been exposed to very hot or very cold temperatures, or
  - been recently swimming or bathing.In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature.
- If ear drops or other ear medications have been placed in the ear canal of one ear, take the temperature in the untreated ear.

## Care and cleaning

The probe tip is the most delicate part of the thermometer. It has to be clean and intact to ensure accurate readings.

If the thermometer is accidentally used without a lens filter, clean the probe tip as follows:

Very gently wipe the surface with a cotton swab or soft cloth moistened with alcohol.

Once the alcohol has dried completely, a new lens filter can be put on and a temperature measurement taken. If the probe tip is damaged, please contact your authorized local service centre.




Use a soft, dry cloth to clean the thermometer display and product exterior. Do not use abrasive cleaners. Never submerge this thermometer in water or any other liquid.

Store the thermometer and lens filters in a dry location free from dust and contamination and away from direct sunlight.

Additional lens filters (LF 40) are available at most stores selling Braun ThermoScan.

For more information, please contact an authorized local Service Centre (see Guarantee Card).

## Trouble-shooting

Situation	Solution
 No lens filter is attached.	Attach new, clean lens filter.
 Measurement before device stabilization	Wait until all the icons stop flashing.
 Ambient temperature is not within the allowed Operating range (50 – 104 °F or 10 – 40 °C).	Allow the thermometer to remain for 30 minutes in a room where the temperature is between 50 and 104 °F or 10 and 40 °C.
 HI = too high  LO = too low Temperature taken is not within typical human temperature range (93.2-108°F or 34-42.2°C)	Make sure the probe tip and lens are clean and a new, clean lens filter is attached. Make sure the thermometer is properly inserted. Then, take a new temperature.
 System error (Er5 - Er9)	Reset the thermometer by removing the battery and then replacing it again, wait 1 minute, then turn on again.  Please contact your authorized Service Center (See Guarantee Card).
 Battery is low.	Insert new battery.
 Battery is too low. Device cannot be powered on to the ready state.	Insert new battery.
Do you have any further questions?	Please contact your authorized Service Center (See Guarantee Card).

## Changing the temperature scale

Your Braun ThermoScan is delivered to you with the Celsius (°C) temperature scale activated. If you wish to switch to Fahrenheit (°F), proceed as follows:

Make sure the thermometer is turned off. Press and hold down the «start» button. After about 8 seconds the display will show «°C» flashing.


Release «start» button and press again and release. The display will show «°F» flashing.


After about 5 seconds the display will auto reset and there will be 2 short beeps.

## Replacing the battery

The thermometer is supplied with one 3V Lithium cell battery (battery type). Insert new battery when the low battery symbol appears on the display.

Open the battery compartment. Remove the battery by inserting the battery cover click into the hole below the battery. Slide battery door back until it snaps in place.

To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

## Calibration

The thermometer is calibrated at the time of manufacture. If at any time you question the accuracy of temperature measurements, please contact an authorized Kaz Service Center.

## Product specifications

Model:	IRT 3020
Displayed temperature range:	34 °C – 42.2 °C (93.2 °F – 108 °F)
Operating ambient temperature range:	10 °C – 40 °C (50 °F – 104 °F)
Display resolution:	0.1 °C or °F
Accuracy for displayed temperature range:	± 0.2 °C (35.5–42 °C) (95.9–107.6 °F) ± 0.3 °C (outside this temperature range)
Clinical repeatability:	± 0.25 °C (± 0.45 °F)
Humidity:	≤ 95% non-condensing
Battery life:	2 years / 1000 measurements

If device is not used within specified temperature and humidity ranges the technical accuracy cannot be ensured.

Type BF equipment

Storage temperature

Operating temperature

See instruction of use

Subject to change without notice.

Standard Reference	Edition	Title
EN 12470-5	2003	Clinical Thermometers – Part 5: Performance of Infrared Ear Thermometers (with maximum device).
EN 60601-1	2006	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.
EN ISO 14971	2007	Medical devices – Application of risk management to medical devices.
EN ISO 10993-1	2009	Biological evaluation of medical devices – Part 1: Evaluation and Testing.
EN 60601-1-2	2007	Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests.

This product conforms to the provisions of the EC directive 93/42/EEC (Medical Device Directive).

Portable and mobile RF communications equipment can affect  
MEDICAL ELECTRICAL EQUIPMENT

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC. For detailed description of EMC requirements please contact an authorized local Service Centre.  
(See Guarantee Card).

Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country.

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