



Braun ThermoScan®
Hygiene caps

SKU: LF40EULA01

- Help to avoid cross-contamination if replaced before each measurement
- Latex and BPA (bisphenol A)-free for improved safety
- Compatible with all ThermoScan[®] Ear thermometers

Features

Fitting a new Braun Hygiene cap each time ensures that dirt or ear wax doesn't build up and clog the sensor window, which could lead to inaccurate temperature readings.

- Help to avoid cross-contamination if replaced before each measurement
- Latex and BPA (bisphenol A)-free for improved safety
- 40 Hygiene caps disposable probe covers per box

In The Box

- 40 Hygiene caps disposable probe covers
- Suitable for all ThermoScan[®] Ear thermometers (IRT6520, IRT6515, IRT6020, IRT4520, IRT4020, IRT3030, IRT3020)



Latex & BPA-free

No irritation

Disposable probe covers are BPA and latex-free, keeping you and your child safe from irritation.



No cross-contamination

Hygienic readings

By replacing the Hygiene cap before each measurement, you avoid the risk of cross-contamination during your readings.